

TOC ANALYSER CALIBRATION WORKSHEETS

(Ref. MICLAB 115)

Inorganic Carbon Standard Worksheet

Date, Name of preparer		line 1
Lot of Na ₂ CO ₃ used in prep.		line 2
Weight of Na ₂ CO ₃ + weigh boat		line 3
Weight of weigh boat		line 4
Weight of Na ₂ CO ₃		line 5

Table 1. **Preparation of the Stock Na₂CO₃ standard**

Concentration of stock solution		line 6
	nominal 500 ppm	

Table 2. **Preparation of the Dilute Na₂CO₃ standard**

Concentration of dilute solution	nominal 25.0 ppm	line 7
----------------------------------	------------------	--------

TOC ANALYSER CALIBRATION WORKSHEETS

(Ref. MICLAB 115)

Organic Carbon Standard Worksheet

Date, Name of preparer		line 1
Lot of KHP used in prep.		line 2
Weight of KHP + weigh boat		line 3
Weight of weigh boat		line 4
Weight of KHP		line 5

Table 1. **Preparation of the Stock KHP standard**

Concentration of stock solution		line 6
	nominal 2000 ppm	

Table 2. **Preparation of the Dilute KHP standard**

Concentration of dilute solution	nominal 25.0 ppm	line 7
----------------------------------	------------------	--------

TOC ANALYSER CALIBRATION WORKSHEETS

(Ref. MICLAB 115)

Calibration Worksheet for Autosampler Calibration

Date, Name of Analyst _____

Measurements of Standards

Average TC of 3rd DI blank		=TC _{DI blank}
Average TC of 3rd TC(KHP) standard		=TC _{ave,TC}
Average TC of 3rd IC (Na ₂ CO ₃) standard		=TC _{ave,IC}
Average IC of 3rd IC (Na ₂ CO ₃) standard		= IC _{ave,IC}
TC value of TC(KHP) standard from formulation		= TC _{std}
Adjusted TC value of the TC(KHP) standard = TC _{std} + TC _{DI blank}		= TC _{adj,TC}

Calculation of Calibration Constants

Current TC calibration constant		=TC _{cal,old}
---------------------------------	--	------------------------

New TC calibration constant	$= \frac{TC_{cal,old} \cdot TC_{adj,TC}}{TC_{ave,TC}}$	=TC _{cal,new}
-----------------------------	--	------------------------

Current IC calibration constant		=IC _{cal,old}
---------------------------------	--	------------------------

New IC calibration constant	$= IC_{cal,old} \cdot \frac{TC_{adj,TC}}{TC_{ave,TC}} \cdot \frac{TC_{ave,IC}}{IC_{ave,IC}}$	=IC _{cal,new}
-----------------------------	--	------------------------

Verification of Calibration

Average TC of 3rd TC (KHP) standard		=TC _{ave,chk,TC}
Average TC of 3rd IC(Na ₂ CO ₃) standard		=TC _{ave,chk,IC}
Average IC of 3rd IC (Na ₂ CO ₃) standard		=IC _{ave,chk,IC}
New TC _{adj,TC} = TC(KHP) _{std} + average TC 3 rd DI Blank =		
TC Std % difference = $\frac{(\text{New TC}_{adj,TC} - TC_{ave,chk,TC})}{\text{New TC}_{adj,TC}} \times 100 =$		(NMT 3%)
IC Std % difference = $(TC_{ave,chk,IC} - IC_{ave,chk,IC}) \times 100 =$		(NMT 3%)

TOC ANALYSER CALIBRATION WORKSHEETS

(Ref. MICLAB 115)

500 ppb Sucrose Standard Worksheet

Date, Name of Preparer

line 1

Lot No. of Sucrose used

line 2

Carbon Assay (% value / 100)

line 3

nominal 0.421

Gross Weight (Sucrose + Boat)

line 4

Tare Weight of Weigh Boat

line 5

Nett Weight of Sucrose

line 6

Carbon Conc. of Stock Solution

line 7

(ppb C = $\frac{\text{g [line 6]} * \text{ca [line 3]} * 10^6}{1 \text{ L}}$)

nominal 25,000 ppb C

Carbon Conc. of Working Solution

line 8

([line 7] ÷ 50)

nominal 500 ppb C

TOC ANALYSER CALIBRATION WORKSHEETS

(Ref. MICLAB 115)

500 ppb Benzoquinone Standard Worksheet

Date, Name of Preparer

line 1

Lot No. of Benzoquinone used

line 2

Carbon Assay (% value / 100)

nominal 0.666

line 3

Gross Weight (Benzoquinone + Boat)

line 4

Tare Weight of Weigh Boat

line 5

Nett Weight of Benzoquinone

nominal 0.075g

line 6

Carbon Conc. of Stock Solution

(ppb C = $\frac{\text{g [line 6]} * \text{ca [line 3]} * 10^6}{1 \text{ L}}$)

nominal 50,000 ppb C

line 7

Carbon Conc. of Working Solution

([line 7] ÷ 100)

nominal 500 ppb C

line 8

TOC ANALYSER CALIBRATION WORKSHEETS

(Ref. MICLAB 115)

Suitability Verification Worksheet

Date and Name of Analyst		line 1
--------------------------	--	--------

Calculated Concentration of Sucrose Standard		line 2
--	--	--------

Calculated Concentration of Suitability Standard		line 3
--	--	--------

Average TOC of Water Blank		line 4
----------------------------	--	--------

Average TOC of Sucrose Standard		line 5
---------------------------------	--	--------

Average TOC of Suitability Standard		line 6
-------------------------------------	--	--------

Response Efficiency $\text{Response Efficiency (\%)} = \frac{(TOC_{\text{suitability}} - TOC_{\text{blank}})}{(TOC_{\text{sucrose}} - TOC_{\text{blank}})} * \frac{C_{\text{sucrose}}}{C_{\text{suitability}}} * 100$	between 85% and 115%	line 7
--	-----------------------------	--------