

# Micromac C

## OPERATING MANUAL



| Version | Prepared: | Verified: | Checked: | Approved: | Issued: |
|---------|-----------|-----------|----------|-----------|---------|
|         | PM        | TM        | QM       | MD        |         |
| 1.3 e   |           |           |          |           |         |

---

### Production Sales and Service.:

Loc. Paduni, 7 – 03012 Anagni (FR) – Tel.: ++39 0775 776058 Fax:++39 0775 772204

mail: [anagni@systea.it](mailto:anagni@systea.it) Sito Web: <http://www.systea.it>



**SYSTEAL Srl**

**$\mu$ MAC- C and  $\mu$ MAC- C MP  
Instrument Maintenance Manual  
Ver. 1.1 e**

**SUMMARY**

|                                                                 |    |
|-----------------------------------------------------------------|----|
| 1. INTRODUCTION .....                                           | 4  |
| 1.1 Analytical cycle.....                                       | 4  |
| 1.2 Calibration cycle .....                                     | 5  |
| 1.3 Technical data .....                                        | 6  |
| 1.4 Analyzer flowchart .....                                    | 11 |
| 1. V6.....                                                      | 11 |
| 2. DESCRIPTION .....                                            | 12 |
| 2.1 LFR reactor .....                                           | 12 |
| 2.2 Sampling .....                                              | 12 |
| 2.3 Analysis, .....                                             | 12 |
| 2.3.1 Reagent injection .....                                   | 13 |
| 2.3.2 Pressure compensation .....                               | 13 |
| 2.3.3 Mixing .....                                              | 13 |
| 2.3.4 Reaction and OD reading.....                              | 14 |
| 2.3.5 Calculations .....                                        | 14 |
| 2.3.6 Wash .....                                                | 14 |
| 2.4 Calibration cycle .....                                     | 15 |
| 3. INSTALLATION .....                                           | 16 |
| 3.1 Physical dimensions and space requirement .....             | 16 |
| 3.2 Dimensional drawing .....                                   | 17 |
| 3.3 Mounting - schematics .....                                 | 18 |
| 3.3.1 Analyzer mounting w/o filter.....                         | 18 |
| 3.3.2 Analyzer mounting with on line self cleaning filter ..... | 19 |
| 3.4 Environmental requirements .....                            | 20 |
| 3.5 Electrical requirements .....                               | 20 |
| 3.6 General check.....                                          | 20 |
| 3.7 Electronic compartment description .....                    | 21 |
| 3.8 Power connection.....                                       | 22 |



|                                                |    |
|------------------------------------------------|----|
| 3.9 Signal input/output.....                   | 22 |
| 4. MAINTENANCE .....                           | 24 |
| 4.1 Reagent and or calibrant replacement ..... | 24 |
| 4.2 Pump tube replacement .....                | 25 |
| 4.3 Cleaning procedure .....                   | 28 |
| 5. TROUBLESHOOTING .....                       | 29 |
| 5.1 Hydraulics .....                           | 29 |
| 5.1.1 Reagents injection.....                  | 29 |
| 6. SPARES LIST .....                           | 33 |



## 1. INTRODUCTION

The **μMAC C** is an automatic chemical analyzer for on-line measurements, based on a new improvement of the exclusive analytical technology developed by SYSTEAL named LFA\* (Loop Flow Analysis);

μMAC C series, allows to run sequentially up to four chemistry on the same analyzer

Possible configurations available for μMAC C series are:

- **MP2:** two chemistries
- **MP3:** three chemistries
- **MP4:** four chemistries

The compactness and the robustness, together with the low reagents consumption, allows the system to be easy used directly for continuous monitoring.

(\*) Patent pending

### 1.1 Analytical cycle

A μMAC C is able to run sequentially one or more analytical cycles (chemistries) under the control of the local microprocessor.

The necessary steps for a complete single analytical cycle are the following:

1. Sample aspiration and introduction in analytical reactor.
2. Sample blank measurement and colorimeter zeroing.
3. Reagents injection with the proper sequence and timing.
4. Sample and reagents mixing.
5. Heating if needed.
6. Optical density increase measurement with a double beam colorimeter.
7. End point measurement storing in non volatile RAM for successive remote and local readings.
8. Concentration value automatic calculation, based on the correlation with internal calibration factor.

The analytical cycle can be automatically repeated for another chemistries, up to four; under the microprocessor control that select a proper parameters set.



## 1.2 **Calibration cycle**

### 1. Calibration with working standard:

The mixed working standard is normally placed inside the  $\mu$ MAC C; it will be used as it is, instead of sample, during calibration cycles. This is the standard procedure to calibrate the analyzer.

### 2. Calibration with stock calibrant and automatic dilution:

An high concentration calibrant is placed inside  $\mu$ MAC C and is connected to one of valves between V1 and V7; to run a calibration cycle, the reactor is filled with distilled water, than a small quantity of the stock calibrant is injected in the reactor (standard addition) and diluted by mixing in the total reactor volume. This procedure is used only for special applications.

For each method, the calibration cycle analyze a standard to calculate the calibration factor; for MP configurations, a mixed calibrant is suggested.



### 1.3 Technical data

Micromac 1000 & Micromac C series - Technical data

| General description   | Data                                                                                                 |
|-----------------------|------------------------------------------------------------------------------------------------------|
| Sensor classification | Colorimetric dual beam with silicon detector or Ion selective electrode (only Micromac 1000 E)       |
| Application           | Seawater, surface water, potable water, waste water, industrial water, soils extract, beverages etc. |

| General specifications                | Data                                                                                                                                                                                                                          |
|---------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Power requirements                    | 12 Vdc, power supply from 220Vac or other local supply to 12 Vdc available as option                                                                                                                                          |
| Humidity                              | Up to 90% not condensable                                                                                                                                                                                                     |
| Temperature range                     | 10 – 40° C Analyzers<br>< 25° Reagents                                                                                                                                                                                        |
| Unit dimensions and weight            | Micromac 1000: 55 x 35 x 11 (d x h x w); 15 kg<br>Micromac 1000 C/E/MP: 80 x 45 x 30 (h x w x d); 25 kg                                                                                                                       |
| Measuring points                      | Micromac 1000, C or E: 2<br>Micromac MP: 1<br>More measuring points available on request                                                                                                                                      |
| Atmospheric pressure range            | No limits                                                                                                                                                                                                                     |
| Effect of electromagnetic fields      | EMC tested according CE compliance                                                                                                                                                                                            |
| Tolerance to electrostatic discharges | EMC tested according CE compliance                                                                                                                                                                                            |
| Pc Specification – O/S                | PC 104 industrial standard Included                                                                                                                                                                                           |
| Positioning and installation details  | Micromac 1000: bench top analyzer, to be installed at around 100 cm from floor, max 2 mt. from sampling point<br>Micromac C/E/MP: wall mounting analyzer to be installed at 60/70 cm from floor, max 2 mt from sampling point |
| LFA reactor volume                    | 10 ml                                                                                                                                                                                                                         |



| General specifications           | Data (cont'd)                                                                                                                                          |
|----------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------|
| Materials in contact with sample | Glass. Silicone, plexiglass, stainless steel AISI 316                                                                                                  |
| CE compliance                    | YES                                                                                                                                                    |
| Year 2000 compliance             | YES                                                                                                                                                    |
| General hazards                  | Only chemical, for details see specific chemistries                                                                                                    |
| Sample conditioning requirements | Depending on the matrix: Imhoff sedimentation stage or special filtering system between 10 and 50 micron in case of presence of not settleable sludge. |

| Sample delivery<br>Operating ranges | Data                                                                                                                                |
|-------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------|
| Temperature range                   | 5° - 30 °C                                                                                                                          |
| Pressure                            | 0.2 – 10 bar; Micromac 1000 max 0.5 bar                                                                                             |
| Flow                                | It depends from the sampling or filtering system; analyzer needs 20 ml of sample for cycle.                                         |
| Connection                          | Micromac C/E/MP: Standard Teflon tube 1.6x3.2 mm other connection on request;<br>Micromac 1000: nipple for silicon tube 1.8 x 4 mm. |
| Waste                               | Micromac C/E/MP: Pressure free 1.8 x 4 mm other on request;<br>Micromac 1000: nipple for silicon tube 1.8 x 4 mm pressure free      |
| Viscosity                           | 5 cps                                                                                                                               |
| Turbidity                           | Not applicable; sample blank correction                                                                                             |
| Colour                              | Not applicable; sample blank correction                                                                                             |
| PH                                  | Depending on the chemical method; normal 3 - 12                                                                                     |



| Signal inputs/outputs     | Data                                                                                                                                                                                                                             |
|---------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 4 – 20 mA – Voltage 0 – 5 | 4 – 20 mA load 400 Ohm, linear response (galvanic isolator module available as option) or 0 – 5 V; separated for each stream and/or parameter in case of Micromac MP                                                             |
| Printer options           | Optional, serial output RS232                                                                                                                                                                                                    |
| Radio or modem links      | Available as option                                                                                                                                                                                                              |
| Earthing details          | Not applicable, 12 Vdc power device                                                                                                                                                                                              |
| Alarms output             | General alarm: 1 potential free switch (SPDT), max load 24V AC DC 0.5 A<br>Limit signal: 1 potential free switch (SPDT), max load 24V AC DC 0.5 A<br>Calibration alarm: 1 potential free switch (SPDT), max load 24V AC DC 0.5 A |
| Serial I/O for signals    | Serial data output RS232 for bidirectional communications, available as option                                                                                                                                                   |
| Analog Input              | Analysis: 1 digital contact with fotocoupler, galvanically isolated<br>Calibration: 1 digital contact with fotocoupler, galvanically isolated                                                                                    |

| Commissioning                                        | Data                                                                                                                                               |
|------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------|
| Manufactures set-up details / pre-installation guide | Available, supplied on order confirmation to allows preparing of installation site, also included in operating manual delivered with the analyzer. |
| Training program                                     | Standard 2 days                                                                                                                                    |
| Factory final test certification                     | Delivered with the analyzer                                                                                                                        |
| Delivery time                                        | 60 days                                                                                                                                            |
| Description of sensor technology                     | Available, enclosed to the Instrument Maintenance Manual                                                                                           |
| Operating sequence                                   | Available, enclosed to the LFA Software Manual                                                                                                     |
| Calibration method                                   | Available, enclosed to the LFA Software Manual                                                                                                     |





| Operational specifications  | Data                                                                                  |
|-----------------------------|---------------------------------------------------------------------------------------|
| Type of analysis            | Batch                                                                                 |
| Analysis interval           | User selectable                                                                       |
| Autodilution                | Available: automatic for off scale samples or manually user selectable                |
| Linearity                   | Over all specified range                                                              |
| Lowest detectable change    | Depending on the method, normal 1% or 2% full scale                                   |
| Repeatability               | Depending on the method, normal better than 1% full scale                             |
| Response time               | From 6 to 30 minutes depending on the method                                          |
| Dead time                   | Not applicable, batch analyzer                                                        |
| Rise time                   | Not applicable, batch analyzer                                                        |
| Fall time                   | Not applicable, batch analyzer                                                        |
| Selectivity / interferences | As for official colorimetric method, or as for Ion selective electrode specifications |
| Drift in 24 hours           | 2% full scale                                                                         |
| Memory effects              | Not present, proper washing procedure                                                 |


| Operational calibration    | Data                                                                                                |
|----------------------------|-----------------------------------------------------------------------------------------------------|
| Frequency / intervals      | Recommended: 24 hours for colorimetric equipment as for manufacturer specifications for electrodes. |
| Single / multi point       | Multipoint: 0 and top of range                                                                      |
| Matrix corrections         | Yes, sample blank correction                                                                        |
| Manual / automatic         | Both                                                                                                |
| Handling / preparation     | None                                                                                                |
| Sample introduction method | Integrated peristaltic pump                                                                         |

| Quality checks                 | Data                                                               |
|--------------------------------|--------------------------------------------------------------------|
| Automatic/manual               | Automatic                                                          |
| Frequency                      | During every calibration cycle, both automatic or manual           |
| Percent of range               | Standard top of range                                              |
| System for control and actions | User selectable deviation warning, from the last calibration value |



| Maintenance                       | Data                                                                                                                                                                                                                           |
|-----------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Frequency schedule                | 2 weeks, reagent replacement if needed, visual inspection                                                                                                                                                                      |
| Expected life of parts            | Pump tube: 1500 hours of running<br>Electrodes: 6 months or as specified from manufacturer<br>Transmission tubing: 1 year                                                                                                      |
| Reagent consumption               | Colorimetric: max 200 microliters for test, usually 100/200 ml of single reagent is enough for one month of operation.<br>Electrodes: buffers or ionic strength adapters as for sample matrix and manufacturer specifications. |
| Reagents stability                | Depending on the method, from 15 to 30 days, temperature < than 25°C; for special applications reagent compartment can be refrigerated on request.                                                                             |
| Manufacturer's service            | Factory or on site service for repairs, special maintenance, upgrading or other; under warranty or not.                                                                                                                        |
| Operator maintenance task         | Reagent replacement, pump tube replacement, cleaning procedures, visual inspections                                                                                                                                            |
| List of consumables               | Pump tube<br>Tubing Kit                                                                                                                                                                                                        |
| Indication of spares availability | 2 weeks delivery time                                                                                                                                                                                                          |

1.4 Analyzer flowchart



## 2. DESCRIPTION

### 2.1 LFR reactor

A simple LFR monoparametric analytical reactor includes several devices (mechanical, hydraulics, optics) sequentially interconnected to form a ring defined LOOP FLOW REACTOR, in short LFR.

A standard LFR hydraulic schematic diagram is shown in the previous drawing; components shown are enough to run up to four chemistries, with the possibility to select two different wavelength. standard hydraulics also allows automatic dilution of off scale samples with a dilution ratio factory selected (usually 10)

Valves V10 S/L allows to open the LFR to pump sample, standard or blank solution in (SAMPLE position) or to keep the LFR closed in LOOP, to allows mixing reaction and so on. Valves Vdil are activated, when a sample dilution is needed.

The cylinder C1 allows the injection of the reagents and the calibrant (only in case of calibration using a stock calibrant see 1.2 Calibration cycle) in the LFR, and contributes to create a strong turbulence inside the circuit, in order to speed up the mixing phase.

The valve V6, when activated, interrupts the LFR hydraulic connection and contributes to the vacuum production inside C1.

The pump P is a single tube peristaltic pump that can be activated in bi-directional direction

The valves VC5 and V8 make the selection of Sample, Diluent and Calibrant.

The cylinder C2 is normally filled with air and contributes to the vacuum production in C1, containing temporarily the volume of liquid pumped out from C1 by the pump P.

### 2.2 Sampling

With valves V10 S/L and VC5 NC, the pump P starts in direct mode; the sample is pumped trough V8, V10, V11, and VC6 in the cylinder C1, and once reached the level LC it goes up on the tube, passes through V6, through the colorimeters flow cells, and than reaches C2 and flows out to the waste through the opposite section of the valves V10 S/L.

After an appropriate sampling time, the cylinder C1 is filled up to the level LH by closing V6 and opening V5.

All the LFR components, including the flow cell, are now filled with sample; the system can read and store the Sample Blank value that will be subtracted from the final OD of the reaction product.

### 2.3 Analysis,

**Following operation are valid for Sample analysis, Reagent Blank and Calibration**

To start the analysis valves V10 S/L are switched NO to close LFR in LOOP.

### **2.3.1 Reagent injection**

To inject reagents into the LFR the valve V6 is closed and the pump P is started in reverse mode until the liquid in C1 will reach the level Ll.

The volume of liquid between LH and Ll has been moved in C2 where has reached the level LV; the strong vacuum produced with this operation in C1 will allows to suck the reagents, with a thin and fast flow-out, activating for a precise time (proportional to the volume), the microvalves V1 or V2 connected to the reagents containers.

### **2.3.2 Pressure compensation**

Opening the valve V6 will restore the original levels in C1 and C2.

To compensate the overpressure, after reagents injection, valves V10 S/L are switched NC for few seconds allowing discharge of excess of liquid contained in the circuit.

### **2.3.3 Mixing**

Then valves V10 S/L are positioned again in LOOP (NO), the pump P is activated in direct mode, to allow a fast mixing of sample and the reagents.

The special position of the IN and OUT tubes inside C1 will contribute to produce a strong turbulence thus allowing a fast mixing and the start of the chemical reaction.



ISO 9001:2008



273

Certificate number AJAEU/10/12116

#### 2.3.4 Reaction and OD reading

The reaction will take place in all the points of the reactor and therefore in the colorimeter flow cell, thus allowing the monitoring of the reaction from time 0 (reagents injection) to the end point.

If methods requiring heating, to speed up the colour development, two types of heating bath can be included in the hydraulics:

- 1) Flow cell heater allowing to heat only the reaction product in the flow cell mainly used for ammonia
- 2) Manifold heater allowing to heat all the reaction product that flows inside the heating bath for a time defined in the analytical cycle.

At the reaction end point the analyzer measures and store the reaction O.D.

#### 2.3.5 Calculations

Reagent Blank OD : stored during calibration and used to calculate the calibration factor

Calibrant OD: stored during calibration and used to calculate the calibration factor

OD Start: Sample blank OD

OD End: Sample OD (Final OD – Sample blank OD)

Sample concentration = OD end x Calibration factor

#### 2.3.6 Wash

Valves V10 S/L are switched NC, all other valves are NO, to allow pumping of diluent (usually distilled water). Pump P starts in direct mode, water fills up LFR washing hydraulics.

After a proper washing time (about 1-2 min) LFR, remains in stand by filled with distilled water ready for a new run.



## 2.4 **Calibration cycle**

Allows the instrument recalibration using a working calibrant corresponding to the higher expected sample value or using a stock calibrant that will be automatically diluted during calibration cycle (see 1.2 calibration cycle)

### 2.4.1.1 **Automatic calibration using working calibrant**

Calibrant is pumped inside the LFR as follows:

with valves V10 S/L and V8 NC, the pump P starts in direct mode; calibrant is pumped through V8, V10, V11 and VC6, in the cylinder C1, and once reached the level LC it goes up on the tube, passes through V6, through the colorimeters flow cells, and then reaches C2 and flows out to the waste through the opposite section of the valves V10 S/L.

After an appropriate sampling time, the cylinder C1 is filled up to the level LH by closing V6 and opening V5.

All the LFR components, including the flow cell, are now filled with sample; the system can read and store the Calibrant Blank value that will be subtracted from the final OD of the reaction product.

Calibrant analysis for calibration cycle follow same operations described in 2.3 Analysis

### 2.4.1.2 **Recalibration using stock calibrants**

The analyzer dilute automatically concentrated calibrant contained in the reagents compartment.

To produce the working calibrant it is used following procedure:

LFR is filled with diluent, valves S/L is switched off and the stock standard is injected using the same procedure used for reagent injection.

The microvalve connected to the calibrant container is activated for a predefined time and, therefore, a precise and reproducible calibrant volume is injected in LFR.

A mixing step starts allowing a fast diffusion of the injected volume in the total LFR volume.

With this technique is possible to obtain a working calibrant up to 200 times less concentrated than the stock calibrant.


The working calibrant just obtained is analysed using the previously described analytical cycle and the O.D. values measured will be used for the calculation of the calibration coefficient.

### 3. INSTALLATION

#### 3.1 Physical dimensions and space requirement


A picture of a  $\mu$ MAC-1000 C MP unit with main subsystems description is shown below; the physical dimensions are the following: 800 (H) x 420 (W) x 275 (D) mm.

Including the reagents and the optional internal battery the  $\mu$ MAC-1000 C MP weights about 25 Kg.





### 3.2 Dimensional drawing




| POSIT. | DESCRIPTION |
|--------|-------------|
| 1      | POWER IN    |
| 2      | SIGNALS OUT |
| 3      | SAMPLE IN   |
| 4      | WASTE       |

| SECT | DESCRIPTION |
|------|-------------|
| A    | TOP         |
| B    | BOTTOM      |

### 3.3 Mounting - schematics


#### 3.3.1 Analyzer mounting w/o filter



**Note: Sample delivery Sample overflos is at customer care**

| Position | Description                                | Position | Description                              |
|----------|--------------------------------------------|----------|------------------------------------------|
| 1        | Sample line                                | 6        | Micromac C                               |
| 2        | Sample overflow                            | 7        | Alarms, signals output and signals input |
| 3        | Overflow waste                             | 8        | 12 V dc line                             |
| 4        | Analyzer waste silicone 2x 4 mm            | 9        | Wall mounting points                     |
| 5        | Sample inlet to analyzer silicone 2 x 4 mm |          |                                          |

### 3.3.2 Analyzer mounting with on line self cleaning filter



| Position | Description                                | Position | Description                                                       |
|----------|--------------------------------------------|----------|-------------------------------------------------------------------|
| 1        | Sample line                                | 7        | Alarms, signals output and signals input                          |
| 2        | On line self Cleaning filter               | 8        | 12 V dc line                                                      |
| 3        | Overflow waste                             | 9        | Wall mounting points                                              |
| 4        | Analyzer waste silicone 2 x 4 mm           | 10       | 12 V Dc line                                                      |
| 5        | Sample inlet to analyzer silicone 2 x 4 mm | 11       | Compressed air                                                    |
| 6        | Micromac C                                 | 12       | Filtered sample pot                                               |
|          |                                            | 13       | Busy contact connection to TB 8 and 9 (see On Line filter manual) |

**Note: Sample delivery to filter inlet is at customer care**



ISO 9001:2008



273

Certificate number AJAEU/10/12116

### **3.4 Environmental requirements**

The system should be operated with the ambient temperature ranging from 10 °C to 40 °C.

A temperature near 0 °C freeze the reagents and the calibrant; a temperature over 40 °C will reduce the reagents stability and life, a Peltier cooled reagents compartment can be provided on request.

### **3.5 Electrical requirements**

μMAC C requires the most safe and commonly diffused 12 VDC power supply; a 220 Vac / 12 Vdc power supply can be provided on request.

The minimum absorption in stand-by is 350 mA, the maximum absorption in operation is 2 A, the mean absorption in operation is 0,7 A.

System is protected against current overload, with a 2 A slow blow fuse placed inside the instrument on the board STA-μMAC-03.

### **3.6 General check**

- Before connecting power supply, open both front doors and inspect for any damage; check if everything is well fixed, board connected, connectors inserted.


### 3.7 Electronic compartment description



**μMAC-1000 C MP Electronic compartment**

The electronic of μMAC-1000 C MP analyser is composed by the following boards:

- **CPU board**, based on a PC-104 architecture, which manages the analyser with a programmable step-by step program, predefined in factory
- **I/O board**, which manages all the LFR actuation and signal acquisition, using the auxiliary RFIO board
- **RFIO board**, the auxiliary input/output board
- **Power supply board**, which gives to the other boards the necessary DC power supplies.




### 3.8 Power connection

Connect the power on the **green connector** of the power board,; take care to connect power cable with right polarity.

For safety reasons it can be useful to measure the current absorption after power on, once connected the AC power; it should be between 350 and 550 mA (if temperature on it can be more than 1 A).

#### POWER SUPPLY BOARD



### 3.9 Signal input/output

#### RFIO board

3.9.1 **Contact port P3A** (left, method 0)


3.9.2 **Contact port P3B** (right, method 1)

3.9.3 **Input contacts P4A positive** (left row)

**Input contacts P4B negative** (right row)



3.9.4 **Output analog P8A** (left, method 0)



**Output analog P8B** (right, method 1)



**NOTE:** Method 0 or Method 1 can be used by the instrument software also to run a single method analyser for Sample Stream 0 and Sample Stream 1 !



| Port | Pin | Internal                                                                                   | Signal type, not act. (actuated)  | External use                                                                     |
|------|-----|--------------------------------------------------------------------------------------------|-----------------------------------|----------------------------------------------------------------------------------|
| P3-A | 1   | Relay K1  | standby(Busy)      Analyser       | To recognise if instrument is busy or if it is ready to accept external commands |
|      | 2   |                                                                                            |                                   |                                                                                  |
|      | 3   |                                                                                            |                                   |                                                                                  |
|      | 4   | Relay K3  | (calibration error)      Method 0 | To recognise if calibration OD was in or out (=error condition) of set limit     |
|      | 5   |                                                                                            |                                   |                                                                                  |
|      | 6   |                                                                                            |                                   |                                                                                  |
|      | 7   | Relay K5  | (calibration error)      Method 1 | To recognise if calibration OD was in or out (=error condition) of set limit     |
|      | 8   |                                                                                            |                                   |                                                                                  |
|      | 9   |                                                                                            |                                   |                                                                                  |
|      | 10  | Relay K7    neg. in                                                                        | standby(Busy)      Method 0       | actuation of optional status                                                     |
|      | S1  | Relay for K7                                                                               | Standby=not actuated              | optional relay for K7                                                            |

|      |    |                                                                                              |                             |                                                                           |
|------|----|----------------------------------------------------------------------------------------------|-----------------------------|---------------------------------------------------------------------------|
| P3-B | 1  | Relay K2   | normal(Error)      Analyser | To recognise if instrument has produced general Error.                    |
|      | 2  |                                                                                              |                             |                                                                           |
|      | 3  |                                                                                              |                             |                                                                           |
|      | 4  | Relay K4  | (High limit)      Method 0  | To recognise if measured concentration is in or out (=alarm) of set limit |
|      | 5  |                                                                                              |                             |                                                                           |
|      | 6  |                                                                                              |                             |                                                                           |
|      | 7  | Relay K6  | (High limit)      Method 1  | To recognise if measured concentration is in or out (=alarm) of set limit |
|      | 8  |                                                                                              |                             |                                                                           |
|      | 9  |                                                                                              |                             |                                                                           |
|      | 10 | Relay K8    neg. in                                                                          | standby(Busy)      Method 1 | actuation of optional status                                              |
|      | S2 | Relay K8                                                                                     | Standby=not actuated        | optional relay for K8                                                     |

|      |     |             |   |                                 |                                                                |
|------|-----|-------------|---|---------------------------------|----------------------------------------------------------------|
| P4-A | 1   | Optokoppler | + | Start analysis (jumper S13B)    | Remote start of analysis cycle, by switching voltage to OC1    |
| P4-B | 1   | 5 - 12 VDC  | - |                                 |                                                                |
| P4-A | 2   | Optokoppler | + | Start calibration (jumper S14B) | Remote start of calibration cycle, by switching voltage to OC2 |
| P4-B | 2   | 5 - 12 VDC  | - |                                 |                                                                |
| P4-A | 3   | Optokoppler | + | not in use                      |                                                                |
| P4-B | 3   | 5 - 12 VDC  | - | (jumpers S15A/S15B/S15C)        |                                                                |
| P4-A | 4   | Optokoppler | + | not in use                      |                                                                |
| P4-B | 4   | 5 - 12 VDC  | - | (jumpers S16A/S16B/S16C)        |                                                                |
| P4-A | 5+6 | GROUND      | - |                                 | Auxiliary supply of external actuation contacts for P4-A/B     |
| P4-B | 5+6 | 12 VDC      | + |                                 |                                                                |

Example: to start the sample from remote, connect P4-B5 with P4-A1 (led of OC1 ON)

|      |   |            |   |                             |                                  |
|------|---|------------|---|-----------------------------|----------------------------------|
| P8-A | 1 | analog out | - | Concentration      Method 0 | 4-20 mA (default)                |
|      | 2 | 4-20 mA    | + | (last measurement)          | alternatively 0-20 mA or 0-5 VDC |
| P8-B | 3 | analog out | - | Concentration      Method 1 | 4-20 mA (default)                |
|      | 4 | 4-20 mA    | + | (last measurement)          | alternatively 0-20 mA or 0-5 VDC |



## 4. MAINTENANCE

### 4.1 Reagent and or calibrant replacement

1. Stop all operations before run any service operation.
2. Check calibrant concentration and measuring unit input.
3. Adit last calibration values (OD Blank and Cal OD) and take note of stored values.
4. Replace reagents and or calibrant solutions, **do not fill old reagents with new solutions**
5. Do not discharge old solutions before completed present procedure
6. Select **Prime** function to fill up reagent lines with new solutions, wait until the analyzer will return in stand by.
7. Select **Calibration** function to run a calibration.
8. Wait until the analyzer will return in stand by.
9. Select **Blank** function to run a reagent blank
10. Wait until the analyzer will return in stand by.
11. Edit the new calibration values (OD Blank and Cal OD) take note of stored values.
12. Start sample analysis

For more information on the reagent replacement procedure and calibration procedure read the software manual on sections:

- Analyzer wash
- Reagents Prime
- Reagent replacement





ISO 9001:2008



273

Certificate number AJAEU/10/12116


## 4.2 Pump tube replacement

**To be performed every 1500 h of pump operation or every six months.**


1. Stop any procedure and run a **Wash** to eliminate any trace of reaction products.
2. Switch the system OFF
3. The pump is placed in the lower front side of the instrument. Unscrew manually the four knurled head screws which hold the pump head. The pumping head is now released and can be opened to remove the silicon pump tube.
4. Replace the tube using the same high density i.e. nipple, keeping in mind that the tube that come out from the lower side of the pump must be connected to the S/L valve, and the tube that come out from the upper side must be connected to the bottom nipple of the cylinder 1.

**Attention: if the tube is installed in the opposite way, there will be severe malfunctioning of the system!**

6. Separate pump halves. Hold pump head as shown, with rollers in the 2, 6, and 10 o'clock positions.



7. Place tubing in the outer port and against the two rollers shown, keeping your thumb on the tubing to hold it in place. Insert tubing key on the back of the rotor shaft and push in as far as possible. Tubing is now positioned deep into the cavity.





ISO 9001:2008




273

Certificate number AJAEU/10/12116

8. With the key firmly pressed against the rotor, turn (counterclockwise) and push down while turning, until tubing has surrounded the rotor.



9. Tubing is now in place. Remove key and position other pump half onto the motor shaft and snap shut. Be careful not to pinch tubing between plastic halves.



10. Remount the pump head, taking firmly joined the two halves, moving eventually the roller block until the shaft aligns with the motor drive.
11. Place the four screws, turning right them alternatively, helping with fingers to maintain them in position and checking the right closing of pump head.
12. Switch on the instrument
13. Run a **Wash** and verify that the diluent flows-in the diluent inlet and flows-out from the waste outlet. If not, switch immediately off the instrument and invert the two connections to the pump tube.  
Restart and repeat the check.
14. If there is no aspiration, there should be some connections not properly done, or the pump tube is not mounted correctly.  
Check the connections for any leakage and eventually repeat the pump tube change procedure.



ISO 9001:2008



273

Certificate number AJAEU/10/12116

**Please note:**

- 1. The pump tubes have a standard functionality of 1500 h with water. The use of reagents needed for the reactions development may decrease the efficiency.**
- 2. The permanence of the system in the STAND BY mode, for some weeks, with no pump movement, may cause the temporary blockage of the pump tube in correspondence of the roll's pressure, blockage which disappear normally after few turns. If the system stop is much more prolonged, in the order of months, it is necessary to remove the tube with the same procedure described for the change, unlock the tube pushing with the fingers. Mount then the tube taking care to the "IN" and "OUT" connections.**
- 3. Pump tube deterioration is not sudden, and is considered as one of the factor determining the recalibration frequency. A tube near to its life limit may need a more frequent recalibration.**



ISO 9001:2008



273

Certificate number AJAEU/10/12116

#### 4.3 Cleaning procedure

1. Switch off the system.
2. Remove, if any, the cadmium reduction column from the  $\mu$ MAC C unit.
3. Place all the lines, sample, diluent and reagents, in the wash solution.
4. Switch on the system.
5. Run Analysis; check that everything is working properly with special attention to the reagents valves injection; the liquid must flow-out very fast and free from the reagents nipples.
6. In case of injection problem refer to troubleshooting chapter.
7. Place all the lines, samples, diluent and reagents, in distilled water.
8. Run analysis to remove any trace of wash solution from the circuit. Repeat twice.
9. Place the new reagents bottles into the reagent box.
10. Start up the analyzer as described in Software Manual.

Please note:

**The cleaning solutions described below are the most commonly used with the  $\mu$ MAC C and must be used when not otherwise specified in the methods description.**

#### CLEANING SOLUTIONS

##### **FOR ACID METHODS:**

|                     |       |        |
|---------------------|-------|--------|
| sodium hypochloride | 7%C12 | 25 ml  |
| sodium hydroxide    |       | 10 gr  |
| distilled water     | q.s.  | 250 ml |

##### **FOR ALKALINE METHODS:**

|                          |      |        |
|--------------------------|------|--------|
| hydrochloric acid, conc. |      | 10 ml  |
| distilled water          | q.s. | 250 ml |



ISO 9001:2008



273

Certificate number AJAEU/10/12116

## 5. TROUBLESHOOTING

### 5.1 Hydraulics

#### 5.1.1 Reagents injection

The reagents and the calibrant are injected in the CYLINDER 1 by in inert valves connected to needles passing through the cylinder cover.

Any liquid injected must flow-out from the needle with a fast and thin flow.

In case of partial obstruction of a reagent tube the liquid will flow-out forming small drops.

If the obstruction is total nothing will be injected.

Before the injection the LFR must produce a vacuum, inside the CYLINDER 1, by transferring a liquid volume from CYLINDER 1 to CYLINDER 2.

The minimum level of liquid transfer to obtain a proper vacuum is 1 cm.



ISO 9001:2008



273

Certificate number AJAEU/10/12116

### 5.1.2 TROUBLESHOOTING GUIDE – TABLE 1

| TROUBLE                            | PROBABLE CAUSE                                                                          | CORRECTIVE ACTION                                                                                                                                                           |
|------------------------------------|-----------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| INJECTION IN DROPS OR NO INJECTION | LIQUID TRANSFER FOR VACUUM PRODUCTION LESS THAN 1 CM                                    | <ul style="list-style-type: none"> <li>- CHECK VALVE V6</li> <li>- CHECK VOLTAGE BELOW 11 VDC</li> <li>- CHANGE PUMP TUBE</li> </ul>                                        |
|                                    | CLOGGING OF THE REAGENT LINE PRESENCE OF PRECIPITATES IN REAGENT CONTAINER              | <ul style="list-style-type: none"> <li>- REMOVE THE REAGENT, PUT THE LINE IN WATER AND START THE WASH PROCEDURE</li> <li>- REPLACE THE REAGENT TUBE</li> </ul>              |
|                                    | SILICON PUMP TUBE BLOCKAGE DUE TO A LONG PERIOD (> 10 DAYS) OF INACTIVITY OF THE SYSTEM | <ul style="list-style-type: none"> <li>- REMOVE THE PUMP TUBE AND CHECK</li> <li>- REPLACE IF NECESSARY</li> </ul>                                                          |
|                                    | REAGENT TUBE PINCHED                                                                    | <ul style="list-style-type: none"> <li>- CHECK THE REAGENT TUBE FROM REAGENT CONTAINER TO THE REAGENT VALVE, CHECK SLEEVINGS</li> <li>- REPLACE THE REAGENT TUBE</li> </ul> |
|                                    | REAGENT VALVE NOT ACTIVATED ELECTRICALLY                                                | <ul style="list-style-type: none"> <li>- CHECK THE ELECTRICAL CONNECTIONS</li> </ul>                                                                                        |
|                                    | REAGENT VALVE PROPERLY ACTIVATED BUT NOT OPENING                                        | <ul style="list-style-type: none"> <li>- CHANGE THE VALVE,</li> <li>- CHECK THE STA-RFIO BOARD AND REPLACE IF NECESSARY</li> </ul>                                          |



ISO 9001:2008



273

Certificate number AJAEU/10/12116

## 5.1.2 TROUBLESHOOTING GUIDE - TABLE 2

| TROUBLE                                             | PROBABLE CAUSE                      | CORRECTIVE ACTION                                                                                                                  |
|-----------------------------------------------------|-------------------------------------|------------------------------------------------------------------------------------------------------------------------------------|
| VALUES AFTER DILUTION NEAR ZERO                     | SILICON TUBE BETWEEN V11 NC BLOCKED | - CHANGE THE TUBE                                                                                                                  |
|                                                     | DILUTION VALVES NOT SWITCHING       | - CHECK THE VALVE WITH AN EXTERNAL 12 VDC POWER SOURCE.<br>- CHECK THE STA-RFIO BOARD AND REPLACE IF NECESSARY                     |
| LEAKAGE IN SOME LFR CONNECTION DUE TO OVER PRESSURE | VOLTAGE SUPPLY OVER 13 VDC          | -CHECK THE EXTERNAL POWER SUPPLY                                                                                                   |
|                                                     | WASTE TUBE PINCHED                  | - CHECK AND FREE THE TUBE                                                                                                          |
|                                                     | VALVE V5 OR V7 NOT ACTIVATED        | - CHECK THE CONNECTIONS AND THE VALVE WITH AN EXTERNAL 12 VDC POWER SOURCE.<br>- CHECK THE STA-RFIO BOARD AND REPLACE IF NECESSARY |
|                                                     | VALVE V5 OR V7 CLOGGED              | - TRY TO CLEAN WITH A SYRINGE CONTAINING A CLEANING SOLUTION<br>- REPLACE VALVES                                                   |



ISO 9001:2008



273

Certificate number AJAEU/10/12116

## 5.1.2 TROUBLESHOOTING GUIDE - TABLE 2 (cont'd)

| TROUBLE                                               | PROBABLE CAUSE                             | CORRECTIVE ACTION                                                                                                                                                                                                                                                                                                                                                                         |
|-------------------------------------------------------|--------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| NO VALUES ON SAMPLE BUT CALIBRATION GIVES GOOD RESULT | VALVE VC5 DEFECTIVE                        | <ul style="list-style-type: none"> <li>- CHECK THAT DURING I SAMPLING THE LIQUID IS ASPIRATED ONLY FROM THE SAMPLE LINE</li> <li>- CHECK THAT DURING SAMPLING NO LIQUID IS ASPIRATED FROM THE WASH AND CALIBRANT LINES</li> <li>- CHECK THE CONNECTIONS AND THE VALVE WITH AN EXTERNAL 12 VDC POWER SOURCE.</li> <li>- CHECK THE STA-RFIO BOARD AND REPLACE IF NECESSARY</li> </ul>       |
| CALIBRATION VALUE IS NOT REPRODUCIBLE                 | VALVE V8 OR VC5 DEFECTIVE                  | <ul style="list-style-type: none"> <li>- CHECK THAT DURING CALIBRANT SAMPLING LIQUID IS ASPIRATED ONLY FORM THE CALIBRATION LINE</li> <li>- CHECK THAT DURING SAMPLING NO LIQUID IS ASPIRATED FROM THE WASH AND SAMPLE LINES</li> <li>- CHECK THE CONNECTIONS AND THE VALVE WITH AN EXTERNAL 12 VDC POWER SOURCE.</li> <li>- CHECK THE STA-RFIO BOARD AND REPLACE IF NECESSARY</li> </ul> |
| BAD VALUES ON SAMPLE AND CALIBRATION                  | ONE OR MORE REAGENTS NOT INJECTED PROPERLY | <ul style="list-style-type: none"> <li>- RUN A REAGENT PRIME AND CHECK THAT ALL VALVES INJECT PROPERLY</li> <li>- RUN THE CLEANING PROCEDURE</li> </ul>                                                                                                                                                                                                                                   |



## 6. SPARES LIST

| P/N          | Description                                                       | Sales unit | Module    | Method                                     | Type      | Suggested stock |
|--------------|-------------------------------------------------------------------|------------|-----------|--------------------------------------------|-----------|-----------------|
| 116-0536-10  | TUBING, T. 0.045                                                  | Mt.        | MANIFOLD  | All                                        | Spare     | 1               |
| 116-0536-15  | TUBING, T. 0.081                                                  | Mt.        | MANIFOLD  | All                                        | Spare     | 1               |
| 116-0537-10  | TUBING, S. 0.045                                                  | Mt.        | MANIFOLD  | All                                        | Spare     | 1               |
| 116-CFLEX-01 | C FLEX TUBING SMALL 1st BUSHING for 562-3019-X1 (6424-50)         | Mt.        | MANIFOLD  | All                                        | Spare     | 1               |
| 116-CFLEX-02 | C FLEX TUBING BIGF 2nd BUSHING 562-3019-x1->178-6275-03 (6424-62) | Mt.        | MANIFOLD  | All                                        | Spare     | 1               |
| 116-R497-110 | TRANSMISSION TUBE, SILICONE 0,045                                 | Mt.        | MANIFOLD  | All                                        | Spare     | 1               |
| 116-R497-150 | TRANSMISSION TUBE, SILICONE 0,065                                 | Mt.        | MANIFOLD  | All                                        | Spare     | 1               |
| 178-6275-03  | REAGENT STRAW TEFLON 1,6X32                                       | PK/10      | MANIFOLD  | All                                        | Spare     | 1               |
| 562-3019-X1  | TUBING TEFLON 1 x 1,5 mm (LFA reag. Tube)                         | Mt.        | MANIFOLD  | All                                        | Spare     | 5               |
| 562-3019-X2  | TUBING TEFLON 1.6 x 3.2 mm (NPA+ Acidic manifold)                 | Mt.        | MANIFOLD  | COD, TP, Total Iron, strong acidic methods | Spare     | 10              |
| 562-3050-X2  | TUBING SILICONE 1.8X4.0MM                                         | Mt.        | MANIFOLD  | All                                        | Spare     | 10              |
| CPU-DISK-00  | DISK ON CHIP FOR CPU-PC104-00                                     | EA/Unità   | HARDWARE  |                                            | On demand |                 |
| CPU-PC104-00 | CPU BOARD WITHOUT DISK ON CHIP                                    | EA/Unità   | HARDWARE  |                                            | On demand |                 |
| MM1-0459-00  | T CONNECTOR - KARTELL                                             | EA/Unità   | MANIFOLD  | All                                        | Spare     | 5               |
| MM1-1250-00  | L CONNECTOR - KARTELL                                             | EA/Unità   | MANIFOLD  | All                                        | Spare     | 5               |
| PRS-SNSR-01  | SAMPLE PRESENCE DETECTOR                                          | EA/Unità   | HARDWARE  | All                                        | On demand |                 |
| QRZ-COIL-02  | QUARTZ COILS UV DIGESTOR 2 mm ID X 4 mm (LFA)                     | EA/Unità   | YDRAULICS | TP, TN; All methods using UV digestor      | Spare     | 1               |

| P/N             | Description                                       | Sales unit | Module    | Method                                                  | Type       | Suggested stock |
|-----------------|---------------------------------------------------|------------|-----------|---------------------------------------------------------|------------|-----------------|
| STA-BATH-01     | HEATING BATH - LFA ANALYZER (TEFLON)              | EA/Unità   | MANIFOLD  | Silica, NO3 Hydrazine, Tot P, Tot N, Tot Fe             | On demand  |                 |
| STA-BATH-15     | HEATING BATH FLOW CELL & 15 MM FLOW CELL ASSY     | EA/Unità   | COLORIM.  | Ammonia Standard range                                  | On demand  |                 |
| STA-BATH-50     | HEATING BATH FLOW CELL & 50 MM FLOW CELL ASSY     | EA/Unità   | COLORIM.  | Ammonia Seawater range                                  | On demand  |                 |
| STA-CDCL-01     | Cd COLUMN LFA ANALYZERS, FILLED                   | EA/Unità   | MANIFOLD  | NO3 Cd reduction                                        | Consumable | 2               |
| STA-DISP-02     | LFA ANALYZER DISPLAY, PARALLEL/GRAPHICS           | EA/Unità   | ARDWARE   |                                                         | On demand  |                 |
| STA-DISPASSY-01 | LFA ANALYZ. DISPLAY ASSY-PAR/GRAPH (KEYP+BD+DISP) | EA/Unità   | HARDWARE  |                                                         | On demand  |                 |
| STA-FER-02      | FERULA SGM TEFLON 1,6x3,2 (Manifold valve)        | EA/Unità   | MANIFOLD  | COD, TP, Strong acidic methods                          | Spare      | 5               |
| STA-FER-03      | FERULA TEFLON 1,0x1,5 (Reagent Valve)             | EA/Unità   | MANIFOLD  | COD, TP, Strong acidic methods                          | Spare      | 5               |
| STA-KEYP-02     | LFA ANALYZER KEYPAD, NEW WITH LOGO                | EA/Unità   | ARDWARE   |                                                         | On demand  |                 |
| STA-M100-00     | REACTION CYLINDER BLOCK                           | EA/Unità   | HYDRAULIC | All                                                     | Spare      | 1               |
| STA-M100-01     | FLOC CELL HOLDER (UNIVERSAL)                      | EA/Unità   | COLORIM.  |                                                         | On demand  |                 |
| STA-M101-01     | FLOW CELL 50 mm                                   | EA/Unità   | COLORIM.  | Seawater applications, High sensitivity methods; TP, TN | Spare      | 1               |
| STA-M101-02     | FLOW CELL 15 mm                                   | EA/Unità   | COLORIM.  | All Standars applications                               | Spare      | 1               |
| STA-M101-09     | FLOW CELL 5 mm                                    | EA/Unità   | COLORIM.  | High level applications                                 | Spare      | 1               |
| STA-M103-01     | PUMP MOTOR                                        | EA/Unità   | HYDRAULIC | All                                                     | On demand  |                 |
| STA-M103-02     | PUMP HEAD                                         | EA/Unità   | HYDRAULIC | All                                                     | Spare      | 1               |

| P/N          | Description                                   | Sales unit | Module    | Method                                       | Type       | Suggested stock |
|--------------|-----------------------------------------------|------------|-----------|----------------------------------------------|------------|-----------------|
| STA-M103-03  | PUMP TUBE, MICROMAC                           | EA/Unità   | MANIFOLD  | All                                          | Consumable | 3               |
| STA-M103-05  | TEFLON PUMP COD REAGENT INJECTION             | EA/Unità   | MANIFOLD  | COD                                          | On demand  |                 |
| STA-M105-05  | 2 WAY VALVE NEW (N.O.)                        | EA/Unità   | MANIFOLD  | All                                          | Spare      | 1               |
| STA-M105-06  | 3 WAY VALVE NEW                               | EA/Unità   | MANIFOLD  | All                                          | Spare      | 2               |
| STA-M105-07  | 2 WAY VALVE NEW (N.C.)                        | EA/Unità   | MANIFOLD  | All                                          | Spare      | 1               |
| STA-M105-12  | 2 WAY VALVE (SMC)                             | EA/Unità   | MANIFOLD  | Reagent valve all analyzers                  |            | 1               |
| STA-M105-14  | 3 WAY VALVE TEFLON (MANIFOLD)                 | EA/Unità   | MANIFOLD  | COD, TP, Strong Acidic methods               | Spare      | 2               |
| STA-M105-15  | 2 WAY VALVE TEFLON (REAGENT INJECTION)        | EA/Unità   | MANIFOLD  | COD, TP, Strong Acidic methods               | Spare      | 2               |
| STA-M200-01  | SAMPLE POT MICROMAC C                         | EA/Unità   | HYDRAULIC | All                                          | On demand  |                 |
| STA-MCOL-00  | COLORIMETER PREAMPLIFIER BOARD                | EA/Unità   | HARDWARE  |                                              | On demand  |                 |
| STA-MPIO-01  | I/O BOARD PC 104                              | EA/Unità   | HARDWARE  |                                              | On demand  |                 |
| STA-NUT-01   | NUT TEFLON 1,6x3,2 - BLACK (Manifold valve)   | EA/Unità   | MANIFOLD  | COD, TP, Strong acidic methods               | Spare      | 5               |
| STA-NUT-03   | NUT TEFLON 1,0x1,5 (Reagent valve)            | EA/Unità   | MANIFOLD  | NPA, COD, TP, Strong acidic methods          | Spare      | 5               |
| STA-ORING-01 | O-RING RILSAN FOR STA-M105-02                 | EA/Unità   | HYDRAULIC |                                              | Spare      | 5               |
| STA-RFIO-04  | DRIVER BOARD NEW (Micromac 1000)              | EA/Unità   | HARDWARE  |                                              |            |                 |
| STA-RFIO-04C | DRIVER BOARD NEW (Micromac C)                 | EA/Unità   | HARDWARE  |                                              | On demand  |                 |
| STA-TSNS-01  | H.B OR F/C TEMPERATURE SENSOR LM35 WITH CABLE | EA/Unità   | COLORIM.  | NH3, SiO2, NO3 Hydrazine, TP, TN Tot Fe etc. | On demand  |                 |

| <b>P/N</b>   | <b>Description</b>                                   | <b>Sales unit</b> | <b>Module</b> | <b>Method</b>                         | <b>Type</b> | <b>Suggested stock</b> |
|--------------|------------------------------------------------------|-------------------|---------------|---------------------------------------|-------------|------------------------|
| STA-uMAC-03E | POWER SUPPLY DRIVER - NEW (graphic/parallel display) | EA/Unità          | HARDWARE      |                                       | On demand   |                        |
| STA-uMAC-06  | KEYBOARD INTERFACE BOARD; PARALL/GRAPH DISPLAY       | EA/Unità          | HARDWARE      |                                       | On demand   |                        |
| STA-UVLM-01  | UV LAMP LFA/FLOWSYS                                  | EA/Unità          | MANIFOLD      | TP, TN; All methods using UV digestor | Consumable  | 2                      |
| STA-UVPW-01  | UV DIGESTOR STARTER                                  | EA/Unità          | HARDWARE      | TP, TN; All methods using UV digestor | Spare       | 1                      |
| STA-WELL-01  | HEATING BATH RESISTOR                                | EA/Unità          | MANIFOLD      |                                       | On demand   |                        |
| STA-XTR110   | ELECTRONIC COMPONENT4- 20 mA                         | EA/Unità          | HARDWARE      |                                       | On demand   |                        |