pH of sample: 3 – 9

0.03 – 1.00 mg/L Cr

Temperature of sample/reagent: 15 – 35°C

Special Notes (For more detailed information: HACH Procedure Manual)

- Please read **Safety Advice** and **Expiration Date** on package.
- **Range of application:** Wastewater, process analysis
- If test is not performed at the **recommended temperature** an **incorrect** result may be obtained.
- **Remark:** The concentration of chromium (III) is obtained mathematically from the difference between chromium (Total) and chromium (VI).

1 – 8

Chromium (Total)

2 , 6 – 8

Chromium (VI)

For Chromium (Total) perform steps **1 – 8**.
For Chromium (VI) perform steps **2** and **6 – 8**.

1

DosiCap™ Zip


Remove foil cap and unscrew DosiCap™ Zip

Carefully remove the foil from the **DosiCap™ Zip** and unscrew cap.

2


Add **2.0 mL** Sample

Pipet **2.0 mL** of sample into the vial.

3


Flip over DosiCap™ Zip and Tighten

Screw the **DosiCap™ Zip** back on.

4


Shake **2 – 3x**

Shake **firmly** 2 – 3 times.

5

Heat **100°C** **1h**


Heat **one hour** at 100°C in the reactor.

6

Apply DosiCap™ B

Cool


Screw an **orange coloured DosiCap™ B** onto the vial.

Note: Do not invert the vial after digestion.
Allow to cool to room temperature.

7

Wait **2 min**

Invert 2 – 3 times.
After **2 minutes** invert again 2 – 3 times.

8

Read

Thoroughly clean the outside of the vial and insert it into the photometer.
The **barcode** is identified, an **automatic evaluation** is carried out after the vial is inserted.

Principle	Interferences		
The digestion converts all chromium present to chromium (VI). Chromium (VI) ions react with 1,5-diphenylcarbazide to form 1,5-diphenylcarbazone, which forms a red complex with chromium (VI).	The ions listed below have been individually checked up to the given concentrations and do not cause interference. We have not determined cumulative effects and the influence of other ions. Larger amounts of iron, copper, and reducing and oxidizing agents give low-bias results. Undissolved chromium is not determined with the determination of chromium(VI). An analyte concentration greatly (above 20 mg/L) in excess of the stated range will adversely affect color formation, resulting in a false reading within the method range. Measurement results can be verified using sample dilutions or standard additions.		
	2000 mg/L: SO ₄ ²⁻ , Na ⁺ , K ⁺ , NO ₃ ⁻	100 mg/L: Mg ²⁺ , NH ₄ ⁺	10 mg/L: Cu ²⁺ , Fe ³⁺
	1000 mg/L: Cl ⁻	50 mg/L: Zn ²⁺ , Ni ²⁺ , Co ²⁺ , Cd ²⁺	5 mg/L: Ag ⁺
	125 mg/L: Ca ²⁺	25 mg/L: Pb ²⁺	1 mg/L: Sn ²⁺

Note: For more detailed information see the HACH Procedure Manual.

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In the U.S.A. – Call toll-free 800-227-4224
Outside the U.S.A. – Contact the HACH office or distributor serving you.