



The art of innovation

Search form

Search

Search

- [About us](#)
- [Products](#)
- [News](#)
- [Contact](#)

Particle Size in accordance with ISO 13322-2:2006

Dynamic particle size analysis in accordance with ISO 13322-2:2006



The Morphious B1 (BT) system is an instrument which measures particles sizes dynamically in accordance with ISO 13322-2:2006 through a moving substrate (belt). For clients interested in ISO 13322-2:2006 particle size analysis, The Morphious B1 (BT) system, is the best choice.

The Morphious B1 (BT) system is an instrument which measures particles sizes dynamically in accordance with ISO 13322-2:2006 through a moving substrate (belt).



Particle size analysis in accordance to ISO 13322-2:2006 are based on a) free falling systems, b) moving substrate (belt) or c) at conveyor discharge point. For cylindrical particles method b) (on moving substrate (belt)) is considered to be the best possible method. The reason for that is that particles are aligned on the belt (meaning in 2 axes) due to their cylindrical shape. Most systems on the market are based on a) (free falling systems). The advantage is that a large number of particles can be measured in a short period of time. Disadvantage is that particles fall and thus are measured in all orientations. Therefore either the particles have to be aligned using a guide plate before passing the camera, or the particles with unwanted orientations have to be identified (by using particle shape parameters) and then excluded from the analysis. For some of samples, guide plates work well – for others they don’t. Then filtering of particles with unwanted orientation becomes necessary – although this seem to work quite well in most cases, the inherent disadvantages of this method are: (i) The number of measured particles is reduced; (b) Particles that are important (e.g. such ones that are much too long) might also be excluded from the analysis. Measuring particles through c) (at conveyer discharge point) does not yield as good results as method b). For clients interested in ISO 13322-2:2006 particle size analysis therefore the The Morphious B1 (BT) system, is the best choice. We can provide in depth information regarding all possible morphometric paramters abstracted from b) moving substrate belt analysis. for more information info@microptik.eu

Technical Specification

Temperature Range	15-40 °C
Power Supply	230V AC/115V AC, 50/60 Hz
Camera	USB 3 CMOS resolution: 1936 x 1216, 161 f/s, Pixel size 5.86 µm: Color calibration module.
Power Consumption	50W
Light Source	LED Ring light, white light spectrum
Size Range	From 1µm till cm
Throughput Resolution (for pellets)	20 µm with 20 mm inspection width, the throughput is 1 – 3 kg/hr 40 µm with 40 mm inspection width, the throughput is 2 – 6 kg/hr 60 µm with 59 mm inspection width, the throughput is 6 – 12 kg/hr
PC	Based on Intel® Core™ i7
Software	Windows®. Powered by Mishell® Mishell is a registered trademark of Microptik BV Windows ia a registered trademark of Microsoft corporation
Dimensions (LxWxH)	50 x 60 x 50 cm
Weight	32 kg

Ask our sales representatives for more information!



[Request for product brochure](#)

- [Particle Size Analysis](#)
 - [Contamination](#)
 - [Size and Shape](#)
 - [Pellets and Powders](#)
 - [Size and Shape V3](#)
 - [Particle Size in accordance with ISO 13322-2:2006](#)
 - [Contamination, Size and Shape](#)
 - [Size and Shape analysis by laser obscuration and video](#)
 - [Liquids](#)
 - [Purity scanner](#)
 - [Gels In Plastic Film](#)
 - [Large Format Film analyzer](#)
 - [Web Inspection Systems](#)
 - [Turn-Key Analyzer Systems](#)

[Heating and Cooling Stages](#)

[Forensic Analysis](#)

[Flash Point Analyzers](#)

[Lab Robot Systems](#)

[3D Laser Topography](#)

[Microscopy](#)

[Consistometers](#)

[XRF Spectrometers](#)

[Spectrometer Systems](#)

[Aerial Exploration](#)

[Software](#)

[Other Products](#)

Log in

Username or e-mail *

Password *

- [Request new password](#)

Log in