



Reference Standards Product Catalog

*International Manufacturer of Reference
Standards and Analytical Reagents*

www.labsertchemical.com

ABOUT LABSERT


Labsert started its operations in order to provide services in laboratory services and technologies field and then company's objective has been turned on to manufacture of certified reference materials and analytical reagents in the course of years. In this direction, infrastructure and quality operations have been completed and firstly, manufacture of analytical inorganic standards and reagents has been initiated. Manufacture of certified reference standard materials which are highly used in spectroscopy and chromatography works, in wet chemistry and internal quality control operations has been implemented. Besides, manufacturing of new reference standards is continued to be included in product portfolio with ongoing research and development activities. Sources allocated for both production quality and research-development activities are increased by annual revisions.

Specialists of basic sciences, engineering, finance and foreign trade as well as professionals having significant experience in quality management, instrumental analysis and laboratory processes take part within the organization of Labsert. The process from manufacture to shipment is carried out in conformity with international quality system and accreditations. Laboratory analysis which is scope of ISO 17025 accreditation is used for certification of reference standards and solutions. We provide our services with a great happiness and pride to national and international universities, research and development centers, analysis laboratories, public institutions as well as many private companies interested in our product portfolio.

Important Differences

- Customer and Distributor Services
- Flexible Manufacturing Process
- Simple "Mix Reference Standard" Search Options
- Product Groups Specific to Country Legislation
- Quality Management
- Price Advantage

QUALITY ASSURANCE



Conformity with International standards

Quality understanding is a matter handled and managed at the highest level within Labsert. This is also a matter used as a management system not only in manufacturing process but also in fields such as management of customer services, shipments, human resources etc.

Labsert has ISO 9001 certification. All analysis and tests during certification works are carried out at ISO 17025 accredited laboratories. Manufacture of the standards is implemented according to ISO Guide 34 conditions.

Analysis Quality of Certification Process

Quality and results of reference standards certification process are among the most important issues for users. As Labsert, we apply a process beyond recognized international certification conditions at this point in order to ensure international quality and to keep reliability of results at the highest level. In addition to generally applied international certification analysis processes:

Raw Materials and Packaging Quality

All raw materials used for manufacture of reference materials are selected at minimum purity level of 98%- 99.999% depending on production field and under possible supply conditions. Attention is paid to ensure that all solvents used for manufacture are ultra-pure (Chromatography grade). Temperature and humidity controls are carried out and recorded daily in manufacturing and storage areas. Low particles LDPE (Low Density Polyethylene) or PPCO (Poly-Propylene Copolymer) bottles are used for product packaging. Following completion of manufacture, bottles are made ready for shipment by being placed in aluminum locked packages to minimize environmental effects and for extra protection.

NIST Traceability

Manufacture of reference standards meets conditions to ensure NIST or SI unit based traceability. During instrumental quantitative analysis within our certification process, calibration curves are generated by NIST traceable certified reference materials. Furthermore, calibration is also carried out by weights having NIST traceability for all mass measurements used for certification.



IMPORTANT DIFFERENCES

- Customer Services

Quality of customer services is a matter handled with the highest importance level within Labsert. Customer services department is not a section established only to make sales, to present offers etc. and far beyond all these, specialists having sufficient experiences and competences to provide consultancy during all stages from the first meeting to end of the service have been selected and assigned. Correct answers to your questions about all theoretical, production, shipment, payment etc. have been presented to you upon discussion with related specialists. You shall be contacted at the latest within 24 hours for all your questions and demands. Thus, you won't have to carry out process-conclusion follow-up of your questions or potential issues and you won't incur loss of time.

- Flexible Manufacturing Process

Manufacturing of Custom-made reference standards as well as stock manufacturing are concluded quickly. You can easily create your special manufacture requests at our web site or consult us about all theoretical matters in relation with your demand.

- Simple "Mix Reference Standard" Search Options

For many users, searching desired reference standard blends at the web site is very difficult as well as being obliged to examine all lists one by one in order to find related mixtures is a matter which causes loss of time. We have included the section "Simple Mix Standard Search" in our web site to avoid it and thus, this difficulty and loss of time has been prevented.

Select Component

Be

Type Name of the component or it's formula and select from the list to add as filter. To delete added component, click on it's name below.

Aluminum (Al) Beryllium (Be)

2 Component as Filter

The suggested concentrations belong to the catalogue product

4 Inorganic Products

Make Custom Request

You can add products listed below, or you can see the details by clicking on it's name or orange button located left side of the products's code.

	Ref.No	Vol.	Product Name	Price
	ICP125.M28003	125ml	125ml-Calibration Standard - 28003	130 €
	ICP125.M28004	125ml	125ml-Calibration Standard - 28004	150 €
	ICP250.M28003	250ml	250ml-Calibration Standard - 28003	182 €
	ICP250.M28004	250ml	250ml-Calibration Standard - 28004	210 €

Description

10 mg/l [Al] Aluminium, 10 mg/l [As] Arsenic, 10 mg/l [Ba] Barium, 10 mg/l [Be] Beryllium, 10 mg/l [Cd] Cadmium, 10 mg/l [Cr] Chromium, 10 mg/l [Co] Cobalt, 10 mg/l [Cu] Copper, 10 mg/l [Pb] Lead, 10 mg/l [Mn] Manganese, 10 mg/l [Ni] Nickel, 10 mg/l [Se] Selenium, 10 mg/l [Ag] Silver, 10 mg/l [Ti] Thallium, 10 mg/l [Th] Thorium, 10 mg/l [U] Uranium, 10 mg/l [V] Vanadium, 10 mg/l [Zn] Zinc

- Product Groups Specific to Country Legislation

Labsert is the first company in the international area to manufacture calibration reference solutions which are specially prepared pursuant to countries' legislative limits and ready to use for analysts. We help users with our calibration solution packages manufactured at 5 concentration points including limit values to establish their analytical device calibration curves without making any preliminary preparation. Furthermore, we may diversify these packages as per customers' special needs.

CONTENTS

AAS Standards & Modifiers

Modifiers, Buffers & Reagents	8
<i>Ionisation Buffers</i>	8
<i>Reagents for AAS Flame</i>	8
<i>Standards for Graphite Furnace</i>	8

ICP & ICP-MS Standards

Single-Element Standards	10
ICP Multi-Element Standards	17
<i>ICP-MS Standards</i>	17 - 24
<i>Instrument Check Standards</i>	17 - 24
<i>Tuning Solutions</i>	17 - 24
Blank Solutions	25

Ion Chromatography Standards

Single-Ion Standards	27
Multi-Ion Standards	29 - 30
Eluent Concentrates	30

CONTENTS

Analytical Reagents

■	pH Standards	32
■	Conductivity Standards	33
■	Turbidity Standards	33
■	TOC Standards	34
■	Miscellaneous Standards	35


Calibration Standards for Legal Regulations

■	Calibration Standards for Water Analysis	37
---	--	----

Cross Reference Standards

■	Cross Reference Standards	39
---	---------------------------	----

AAS STANDARDS



AAS Standards and Reagents

AAS Standards and Reagents

Ionisation Buffers

Product Code	Symbol	Chapter	Matrix	Volume	Value
AAS125.SR20000.2NA	Sr	AAS 2% Strontium Ionisation Buffer Standard in 2% HNO ₃	%2 HNO ₃	125 ml	2% Sr
AAS125.CS20000.1NA	Cs	AAS 2% Cesium Ionisation Buffer Standard in 1% HNO ₃	%1 HNO ₃	125 ml	2% Cs
AAS125.LI20000.1NA	Li	AAS 2% Lithium Ionisation Buffer Standard in 1% HNO ₃	%1 HNO ₃	125 ml	2% Li

Reagents for AAS Flame

Product Code	Symbol	Chapter	Matrix	Volume	Value
AAS125.LA20000.2CA	La	AAS 2% Lanthanum Reagents for AAS Flame in 2% HCl	%2 HCl	125 ml	2% La
AAS500.LA20000.2CA	La	AAS 2% Lanthanum Reagents for AAS Flame in 2% HCl	%2 HCl	500 ml	2% La
AAS125.SR20000.2CA	Sr	AAS 2% Strontium Reagents for AAS Flame in 2% HCl	%2 HCl	125 ml	2% Sr
AAS500.SR20000.2CA	Sr	AAS 2% Strontium Reagents for AAS Flame in 2% HCl	%2 HCl	500 ml	2% Sr

Reference Standards for Graphite Furnace

Product Code	Symbol	Chapter	Matrix	Volume	Value
AAS125.SB002.1NA	Sb	AAS 0.02 mg/l Antimony Calibration Standard in 1% HNO ₃	%1 HNO ₃	125 ml	0.02 mg/l
AAS125.AR002.1NA	As	AAS 0.02 mg/l Arsenic Calibration Standard in 1% HNO ₃	%1 HNO ₃	125 ml	0.02 mg/l
AAS125.CD002.1NA	Cd	AAS 0.02 mg/l Cadmium Calibration Standard in 1% HNO ₃	%1 HNO ₃	125 ml	0.02 mg/l
AAS125.PB002.1NA	Pb	AAS 0.02 mg/l Lead Calibration Standard in 1% HNO ₃	%1 HNO ₃	125 ml	0.02 mg/l
AAS125.MG002.1NA	Mg	AAS 0.02 mg/l Magnesium Calibration Standard in 1% HNO ₃	%1 HNO ₃	125 ml	0.02 mg/l
AAS125.NI002.1NA	Ni	AAS 0.02 mg/l Nickel Calibration Standard in 1% HNO ₃	%1 HNO ₃	125 ml	0.02 mg/l

ICP-MS STANDARDS



ICP-MS Single Standards

ICP-MS Single Standards

Product Code	Symbol	Description	Matrix	Volume	Value
ICP125.SR100.2NA	Sr	ICP-MS 100 mg/L Strontium standard in %2 HNO3	%2 HNO3	125 ml	100 mg/L
ICP125.V100.2NA	V	ICP-MS 100 mg/L Vanadium standard in %2 HNO3	%2 HNO3	125 ml	100 mg/L
ICP250.SR1000.2NA	Sr	ICP-MS 1000 mg/L Strontium Standard in %2 HNO3	%2 HNO3	250 ml	1000 mg/L
ICP500.SR1000.2NA	Sr	ICP-MS 1000 mg/L Strontium Standard in %2 HNO3	%2 HNO3	500 ml	1000 mg/L
ICP250.V1000.2NA	V	ICP-MS 1000 mg/L Vanadium Standard in %2 HNO3	%2 HNO3	250 ml	1000 mg/L
ICP500.V1000.2NA	V	ICP-MS 1000 mg/L Vanadium Standard in %2 HNO3	%2 HNO3	500 ml	1000 mg/L
ICP125.NA100.01NA	Na	ICP-MS 100 mg/L Sodium Standard in %0.1 HNO3	%0.1 HNO3	125 ml	100 mg/L
ICP250.NA1000.2NA	Na	ICP-MS 1000 mg/L Sodium Standard in %0.1 HNO3	%0.1 HNO3	250 ml	1000 mg/L
ICP500.NA1000.01NA	Na	ICP-MS 1000 mg/L Sodium Standard in %0.1 HNO3	%0.1 HNO3	500 ml	1000 mg/L
ICP125.AG100.5NA	Ag	ICP-MS 100 mg/L Silver Standard in %5 HNO3	%5 HNO3	125 ml	100 mg/L
ICP250.AG1000.5NA	Ag	ICP-MS 1000 mg/L Silver Standard in %5 HNO3	%5 HNO3	250 ml	1000 mg/L
ICP500.AG1000.5NA	Ag	ICP-MS 1000 mg/L Silver Standard in %5 HNO3	%5 HNO3	500 ml	1000 mg/L
ICP125.K100.01NA	K	ICP-MS 100 mg/L Potassium Standard in %0.1 HNO3	%0.1 HNO3	125 ml	100 mg/L
ICP250.K1000.01NA	K	ICP-MS 1000 mg/L Potassium Standard in %0.1 HNO3	%0.1 HNO3	250 ml	1000 mg/L
ICP500.K1000.01NA	K	ICP-MS 1000 mg/L Potassium Standard in %0.1 HNO3	%0.1 HNO3	500 ml	1000 mg/L
ICP125.Y100.2NA	Y	ICP-MS 100 mg/L Yttrium Standard in %2 HNO3	%2 HNO3	125 ml	100 mg/L
ICP250.Y1000.2NA	Y	ICP-MS 1000 mg/L Yttrium standard in %2 HNO3	%2 HNO3	250 ml	1000 mg/L
ICP500.Y1000.2NA	Y	ICP-MS 1000 mg/L Yttrium Standard in %2 HNO3	%2 HNO3	500 ml	1000 mg/L
ICP125.MG100.01NA	Mg	ICP-MS 100 mg/L Magnesium Standard in %0.1 HNO3	%0.1 HNO3	125 ml	100 mg/L
ICP250.MG1000.01NA	Mg	ICP-MS 1000 mg/L Magnesium Standard in %0.1 HNO3	%0.1 HNO3	250 ml	1000 mg/L
ICP500.MG1000.3NA	Mg	ICP-MS 1000 mg/L Magnesium Standard in %0.1 HNO3	%0.1 HNO3	500 ml	1000 mg/L
ICP125.LI100.01NA	Li	ICP-MS 100 mg/L Lithium Standard in %0.1 HNO3	%0.1 HNO3	125 ml	100 mg/L
ICP250.LI1000.01NA	Li	ICP-MS 1000 mg/L Lithium Standard in %0.1 HNO3	%0.1 HNO3	250 ml	1000 mg/L
ICP500.LI1000.01NA	Li	ICP-MS 1000 mg/L Lithium Standard in %0.1 HNO3	%0.1 HNO3	500 ml	1000 mg/L
ICP125.CE100.3NA	Ce	ICP-MS 100 mg/L Cerium Standard in %3 HNO3	%3 HNO3	125 ml	100 mg/L
ICP250.CE1000.3NA	Ce	ICP-MS 1000 mg/L Cerium Standard in %3 HNO3	%3 HNO3	250 ml	1000 mg/L
ICP500.CE1000.3NA	Ce	ICP-MS 1000 mg/L Cerium Standard in %3 HNO3	%3 HNO3	500 ml	1000 mg/L
ICP125.NI100.2NA	Ni	ICP-MS 100 mg/L Nickel Standard in %2 HNO3	%2 HNO3	125 ml	100 mg/L
ICP250.NI1000.2NA	Ni	ICP-MS 1000 mg/L Nickel Standard in %2 HNO3	%2 HNO3	250 ml	1000 mg/L
ICP500.NI1000.2NA	Ni	ICP-MS 1000 mg/L Nickel Standard in %2 HNO3	%2 HNO3	500 ml	1000 mg/L
ICP125.HG100.2CA	Hg	ICP-MS 100 mg/L Mercury Standard in %2 HCl	%2 HCl	125 ml	100 mg/L
ICP250.HG1000.2CA	Hg	ICP-MS 1000 mg/L Mercury Standard in %2 HCl	%2 HCl	250 ml	1000 mg/L
ICP500.HG1000.2CA	Hg	ICP-MS 1000 mg/L Mercury Standard in %2 HCl	%2 HCl	500 ml	1000 mg/L
ICP125.PB100.05NA	Pb	ICP-MS 100 mg/L Lead in %0.5 HNO3	%0.5 HNO3	125 ml	100 mg/L
ICP250.PB1000.05NA	Pb	ICP-MS 1000 mg/L Lead Standard in %0.5 HNO3	%0.5 HNO3	250 ml	1000 mg/L
ICP500.PB1000.05NA	Pb	ICP-MS 1000 mg/L Lead Standard in %0.5 HNO3	%0.5 HNO3	500 ml	1000 mg/L
ICP125.FE100.3NA	Fe	ICP-MS 100 mg/L Iron Standard in %3 HNO3	%3 HNO3	125 ml	100 mg/L
ICP250.FE1000.3NA	Fe	ICP-MS 1000 mg/L Iron Standard in %3 HNO3	%3 HNO3	250 ml	1000 mg/L
ICP500.FE1000.3NA	Fe	ICP-MS 1000 mg/L Iron Standard in %3 HNO3	%3 HNO3	500 ml	1000 mg/L

ICP-MS Single Standards

ICP-MS Single Standards

ICP125.CO100.2NA	Co	ICP-MS 100 mg/L Cobalt Standard in %2 HNO3	%2 HNO3	125 ml	100 mg/L
ICP250.CO1000.2NA	Co	ICP-MS 1000 mg/L Cobalt Standard in %2 HNO3	%2 HNO3	250 ml	1000 mg/L
ICP500.CO1000.2NA	Co	ICP-MS 1000 mg/L Cobalt Standard in %2 HNO3	%2 HNO3	500 ml	1000 mg/L
ICP125.CR100.2NA	Cr	ICP-MS 100 mg/L Chromium Standard in %2 HNO3	%2 HNO3	125 ml	100 mg/L
ICP250.CR1000.2NA	Cr	ICP-MS 1000 mg/L Chromium Standard in %2 HNO3	%2 HNO3	250 ml	1000 mg/L
ICP500.CR1000.2NA	Cr	ICP-MS 1000 mg/L Chromium Standard in %2 HNO3	%2 HNO3	500 ml	1000 mg/L
ICP125.CD100.3NA	Cd	ICP-MS 100 mg/L Cadmium in %3 HNO3	%3 HNO3	125 ml	100 mg/L
ICP250.CD1000.3NA	Cd	ICP-MS 1000 mg/L Cadmium Standard in %3 HNO3	%3 HNO3	250 ml	1000 mg/L
ICP500.CD1000.3NA	Cd	ICP-MS 1000 mg/L Cadmium Standard in %3 HNO3	%3 HNO3	500 ml	1000 mg/L
ICP125.B100.01NA	B	ICP-MS 100 mg/L Boron Standard in %0.1 HNO3	%0.1 HNO3	125 ml	100 mg/L
ICP250.B1000.01NA	B	ICP-MS 1000 mg/L Boron Standard in %0.1 HNO3	%0.1 HNO3	250 ml	1000 mg/L
ICP500.B1000.01NA	B	ICP-MS 1000 mg/L Boron Standard in %0.1 HNO3	%0.1 HNO3	500 ml	1000 mg/L
ICP125.AS100.2CA	As	ICP-MS 100 mg/L Arsenic Standard in %2 HCl	%2 HCl	125 ml	100 mg/L
ICP250.AS1000.2CA	As	ICP-MS 1000 mg/L Arsenic Standard in %2 HCl	%2 HCl	250 ml	1000 mg/L
ICP500.AS1000.2CA	As	ICP-MS 1000 mg/L Arsenic Standard in %2 HCl	%2 HCl	500 ml	1000 mg/L
ICP125.AL100.1NA	Al	ICP-MS 100 mg/L Aluminium Standard in %1 HNO3	%1 HNO3	125 ml	100 mg/L
ICP250.AL1000.1NA	Al	ICP-MS 1000 mg/L Aluminum Standard in %1 HNO3	%1 HNO3	250 ml	1000 mg/L
ICP500.AL1000.1NA	Al	ICP-MS 1000 mg/L Aluminum Standard in %1 HNO3	%1 HNO3	500 ml	1000 mg/L
ICP125.TB100.3NA	Tb	ICP-MS 100 mg/L Terbium Standard in %3 HNO3	%3 HNO3	125 ml	100 mg/L
ICP250.TB1000.3NA	Tb	ICP-MS 1000 mg/L Terbium Standard in %3 HNO3	%3 HNO3	250 ml	1000 mg/L
ICP500.TB1000.3NA	Tb	ICP-MS 1000 mg/L Terbium Standard in %3 HNO3	%3 HNO3	500 ml	1000 mg/L
ICP125.SC100.2NA	Sc	ICP-MS 100 mg/L Scandium Standard in %2 HNO3	%2 HNO3	125 ml	100 mg/L
ICP125.SC1000.2NA	Sc	ICP-MS1000 mg/L Scandium Standard in %2 HNO3	%2 HNO3	250 ml	1000 mg/L
ICP125.RH100.3NA	Rh	ICP-MS 100 mg/L Rhodium Standard in %3 HNO3	%3 HNO3	125 ml	100 mg/L
ICP125.RH1000.3NA	Rh	ICP-MS 1000 mg/L Rhodium Standard in %3 HNO3	%3 HNO3	250 ml	1000 mg/L
ICP125.LU100.3NA	Lu	ICP-MS 100 mg/L Lutetium Standard in %3 HNO3	%3 HNO3	125 ml	100 mg/L
ICP250.LU1000.3NA	Lu	ICP-MS 1000 mg/L Lutetium Standard in %3 HNO3	%3 HNO3	250 ml	1000 mg/L
ICP500.LU1000.3NA	Lu	ICP-MS 1000 mg/L Lutetium Standard in %3 HNO3	%3 HNO3	500 ml	1000 mg/L
ICP125.GE100.1NA	Ge	ICP-MS 100 mg/L Germanium Standard in %1 HNO3	%1 HNO3	125 ml	100 mg/L
ICP125.GE1000.1NA	Ge	ICP-MS 1000 mg/L Germanium standard in %1 HNO3	%1 HNO3	250 ml	1000 mg/L
ICP125.MN100.3NA	Mn	ICP-MS 100 mg/L Manganese Standard in %3 HNO3	%3 HNO3	125 ml	100 mg/L
ICP250.MN1000.3NA	Mn	ICP-MS 1000 mg/L Manganese Standard in %3 HNO3	%3 HNO3	250 ml	1000 mg/L
ICP500.MN1000.3NA	Mn	ICP-MS 1000 mg/L Manganese Standard in %3 HNO3	%3 HNO3	500 ml	1000 mg/L
ICP125.SE100.2NA	Se	ICP-MS 100 mg/L Selenium Standard in %2 HNO3	%2 HNO3	125 ml	100 mg/L
ICP250.SE1000.2NA	Se	ICP-MS 1000 mg/L Selenium Standard in %2 HNO3	%2 HNO3	250 ml	1000 mg/L
ICP500.SE1000.2NA	Se	ICP-MS 1000 mg/L Selenium Standard in %2HNO3	%2HNO3	500 ml	1000 mg/L
ICP125.BE100.3NA	Be	ICP-MS 100 mg/L Beryllium Standard in %3 HNO3	%3 HNO3	125 ml	100 mg/L
ICP250.BE1000.3NA	Be	ICP-MS 1000 mg/L Beryllium Standard in %3 HNO3	%3 HNO3	250 ml	1000 mg/L
ICP500.BE1000.3NA	Be	ICP-MS 1000 mg/L Beryllium Standard in %3 HNO3	%3 HNO3	500 ml	1000 mg/L
ICP125.SI100.3NA	Si	ICP-MS 100mg/L Silicon Standard in %3 HNO3	%3 HNO3	125 ml	100 mg/L

ICP-MS Single Standards

ICP-MS Single Standards

ICP250.SI1000.3NA	Si	ICP-MS 1000 mg/L Silicon Standard in %3 HNO3	%3 HNO3	250 ml	1000 mg/L
ICP500.SI1000.3NA	Si	ICP-MS 1000 mg/L Silicon Standard in %3 HNO3	%3 HNO3	500 ml	1000 mg/L
ICP125.TI100.3NA	Ti	ICP-MS 100mg/L Titanium Standard in %3 HNO3	%3 HNO3	125 ml	100 mg/L
ICP250.TI1000.3NA	Ti	ICP-MS 1000 mg/L Titanium Standard in %3 HNO3	%3 HNO3	250 ml	1000 mg/L
ICP500.TI1000.3NA	Ti	ICP-MS 1000 mg/L Titanium Standard in %3 HNO3	%3 HNO3	500 ml	1000 mg/L
ICP125.BA100.2NA	Ba	ICP-MS 100 mg/L Barium Standard in %2 HNO3	%2 HNO3	125 ml	100 mg/L
ICP250.BA1000.2NA	Ba	ICP-MS 1000 mg/L Barium Standard in %2 HNO3	%2 HNO3	250 ml	1000 mg/L
ICP500.BA1000.2NA	Ba	ICP-MS 1000 mg/L Barium Standard in %2 HNO3	%2 HNO3	500 ml	1000 mg/L
ICP125.ZN100.2NA	Zn	ICP-MS 100mg/L Zinc Standard in %2 HNO3	%2 HNO3	125 ml	100 mg/L
ICP250.ZN1000.2NA	Zn	ICP-MS 1000 mg/L Zinc Standard in %2 HNO3	%2 HNO3	250 ml	1000 mg/L
ICP500.ZN1000.2NA	Zn	ICP-MS 1000 mg/L Zinc Standard in %2 HNO3	%2 HNO3	500 ml	1000 mg/L
ICP125.CU100.2NA	Cu	ICP-MS 100mg/L Copper Standard in %2 HNO3	%2 HNO3	125 ml	100 mg/L
ICP250.CU1000.2NA	Cu	ICP-MS 1000 mg/L Copper Standard in %2 HNO3	%2 HNO3	250 ml	1000 mg/L
ICP500.CU1000.2NA	Cu	ICP-MS 1000 mg/L Copper Standard in %2 HNO3	%2 HNO3	500 ml	1000 mg/L
ICP125.TI100.2NA	Tl	ICP-MS 100mg/L Thallium Standard in %2 HNO3	%2 HNO3	125 ml	100 mg/L
ICP250.TI1000.2NA	Tl	ICP-MS 1000 mg/L Thallium Standard in %2 HNO3	%2 HNO3	250 ml	1000 mg/L
ICP500.TI1000.2NA	Tl	ICP-MS 1000 mg/L Thallium Standard in %2 HNO3	%2 HNO3	500 ml	1000 mg/L
ICP125.IN100.3NA	In	ICP-MS 100mg/L Indium Standard in %3 HNO3	%3 HNO3	125 ml	100 mg/L
ICP250.IN1000.3NA	In	ICP-MS 1000 mg/L Indium Standard in %3 HNO3	%3 HNO3	250 ml	1000 mg/L
ICP500.IN1000.3NA	In	ICP-MS 1000 mg/L Indium Standard in %3 HNO3	%3 HNO3	500 ml	1000 mg/L
ICP125.CA100.01NA	Ca	ICP-MS 100mg/L Calcium Standard in %2 HNO3	%2 HNO3	125 ml	100 mg/L
ICP250.CA1000.01NA	Ca	ICP-MS 1000 mg/L Calcium Standard in %2 HNO3	%2 HNO3	250 ml	1000 mg/L
ICP500.CA1000.01NA	Ca	ICP-MS 1000 mg/L Calcium Standard in %2 HNO3	%2 HNO3	500 ml	1000 mg/L
ICP125.BI100.3NA	Bi	ICP-MS 100 mg/L Bismuth Standard in %3 HNO3	%3 HNO3	125 ml	100 mg/L
ICP250.BI1000.3NA	Bi	ICP-MS 1000 mg/L Bismuth Standard in %3 HNO3	%3 HNO3	250 ml	1000 mg/L
ICP500.BI1000.3NA	Bi	ICP-MS 1000 mg/L Bismuth Standard in %3 HNO3	%3 HNO3	500 ml	1000 mg/L
ICP125.CS100.2NA	Cs	ICP-MS 100 mg/L Cesium Standard in %2 HNO3	%2 HNO3	125 ml	100 mg/L
ICP250.CS1000.2NA	Cs	ICP-MS 1000 mg/L Cesium Standard in %2 HNO3	%2 HNO3	250 ml	1000 mg/L
ICP500.CS1000.2NA	Cs	ICP-MS 1000 mg/L Cesium Standard in %2 HNO3	%2 HNO3	500 ml	1000 mg/L
ICP125.AU100.2CA	Au	ICP-MS 100 mg/L Gold Standard in %2 HCl	%2 HCl	125 ml	100 mg/L
ICP250.AU1000.2CA	Au	ICP-MS 1000 mg/L Gold Standard in %2 HCl	%2 HCl	250 ml	1000 mg/L
ICP500.AU1000.2CA	Au	ICP-MS 1000 mg/L Gold Standard in %2 HCl	%2 HCl	500 ml	1000 mg/L
ICP125.PD100.2NA	Pd	ICP-MS 100 mg/L Palladium Standard in %2 HNO3	%2 HNO3	125 ml	100 mg/L
ICP250.PD1000.2NA	Pd	ICP-MS 1000 mg/L Palladium Standard in %2 HNO3	%2 HNO3	250 ml	1000 mg/L
ICP500.PD1000.2NA	Pd	ICP-MS 1000 mg/L Palladium Standard in %2 HNO3	%2 HNO3	500 ml	1000 mg/L
ICP125.PT100.2CA	Pt	ICP-MS 100 mg/L Platinum Standard in %2 HCl	%2 HCl	125 ml	100 mg/L
ICP250.PT1000.2CA	Pt	ICP-MS 1000 mg/L Platinum Standard in %2 HCl	%2 HCl	250 ml	1000 mg/L
ICP500.PT1000.2CA	Pt	ICP-MS 1000 mg/L Platinum Standard in %2 HCl	%2 HCl	500 ml	1000 mg/L
ICP125.RU100.2CA	Ru	ICP-MS 100 mg/L Ruthenium Standard in %2 HCl	%2 HCl	125 ml	100 mg/L
ICP250.RU1000.2CA	Ru	ICP-MS 1000 mg/L Ruthenium Standard in %2 HCl	%2 HCl	250 ml	1000 mg/L

ICP-MS Single Standards

ICP-MS Single Standards

ICP500.RU1000.2CA	Ru	ICP-MS 1000 mg/L Ruthenium Standard in %2 HCl	%2 HCl	500 ml	1000 mg/L
ICP125.RB100.2NA	Rb	ICP-MS 100 mg/L Rubidium Standard in %2 HNO3	%2 HNO3	125 ml	100 mg/L
ICP250.RB1000.2NA	Rb	ICP-MS 1000 mg/L Rubidium Standard in %2 HNO3	%2 HNO3	250 ml	1000 mg/L
ICP500.RB1000.2NA	Rb	ICP-MS 1000 mg/L Rubidium Standard in %2 HNO3	%2 HNO3	500 ml	1000 mg/L
ICP125.U100.2NA	U	ICP-MS 100 mg/L Uranium Standard in %2 HNO3	%2 HNO3	125 ml	100 mg/L
ICP250.U1000.2NA	U	ICP-MS 1000 mg/L Uranium Standard in %2 HNO3	%2 HNO3	250 ml	1000 mg/L
ICP500.U1000.2NA	U	ICP-MS 1000 mg/L Uranium Standard in %2 HNO3	%2 HNO3	500 ml	1000 mg/L
ICP125.DY100.2NA	Dy	ICP-MS 100 mg/L Dysprosium Standard in %2 HNO3	%2 HNO3	125 ml	100 mg/L
ICP250.DY1000.2NA	Dy	ICP-MS 1000 mg/L Dysprosium Standard in %2 HNO3	%2 HNO3	250 ml	1000 mg/L
ICP500.DY1000.2NA	Dy	ICP-MS 1000 mg/L Dysprosium Standard in %2 HNO3	%2 HNO3	500 ml	1000 mg/L
ICP125.ER100.2NA	Er	ICP-MS 100 mg/L Erbium Standard in %2 HNO3	%2 HNO3	125 ml	100 mg/L
ICP250.ER1000.2NA	Er	ICP-MS 1000 mg/L Erbium Standard in %2 HNO3	%2 HNO3	250 ml	1000 mg/L
ICP500.ER1000.2NA	Er	ICP-MS 1000 mg/L Erbium Standard in %2 HNO3	%2 HNO3	500 ml	1000 mg/L
ICP125.EU100.2NA	Eu	ICP-MS 100 mg/L Europium Standard in %2 HNO3	%2 HNO3	125 ml	100 mg/L
ICP250.EU1000.2NA	Eu	ICP-MS 1000 mg/L Europium Standard in %2 HNO3	%2 HNO3	250 ml	1000 mg/L
ICP500.EU1000.2NA	Eu	ICP-MS 1000 mg/L Europium Standard in %2 HNO3	%2 HNO3	500 ml	1000 mg/L
ICP125.LA100.2NA	La	ICP-MS 100 mg/L Lanthanum Standard in %2 HNO3	%2 HNO3	125 ml	100 mg/L
ICP250.LA1000.2NA	La	ICP-MS 1000 mg/L Lanthanum Standard in %2 HNO3	%2 HNO3	250 ml	1000 mg/L
ICP500.LA1000.2NA	La	ICP-MS 1000 mg/L Lanthanum Standard in %2 HNO3	%2 HNO3	500 ml	1000 mg/L
ICP125.PR100.2NA	Pr	ICP-MS 100 mg/L Praseodymium Standard in %2 HNO3	%2 HNO3	125 ml	100 mg/L
ICP250.PR1000.2NA	Pr	ICP-MS 1000 mg/L Praseodymium Standard in %2 HNO3	%2 HNO3	250 ml	1000 mg/L
ICP500.PR1000.2NA	Pr	ICP-MS 1000 mg/L Praseodymium Standard in %2 HNO3	%2 HNO3	500 ml	1000 mg/L
ICP125.GD100.2NA	Gd	ICP-MS 100 mg/L Gadolinium Standard in %2 HNO3	%2 HNO3	125 ml	100 mg/L
ICP250.GD1000.2NA	Gd	ICP-MS 1000 mg/L Gadolinium Standard in %2 HNO3	%2 HNO3	250 ml	1000 mg/L
ICP500.GD1000.2NA	Gd	ICP-MS 1000 mg/L Gadolinium Standard in %2 HNO3	%2 HNO3	500 ml	1000 mg/L
ICP125.HO100.2NA	Ho	ICP-MS 100 mg/L Holmium Standard in %2 HNO3	%2 HNO3	125 ml	100 mg/L
ICP250.HO1000.2NA	Ho	ICP-MS 1000 mg/L Holmium Standard in %2 HNO3	%2 HNO3	250 ml	1000 mg/L
ICP500.HO1000.2NA	Ho	ICP-MS 1000 mg/L Holmium Standard in %2 HNO3	%2 HNO3	500 ml	1000 mg/L
ICP125.TM100.2NA	Tm	ICP-MS 100 mg/L Thulium Standard in %2 HNO3	%2 HNO3	125 ml	100 mg/L
ICP250.TM1000.2NA	Tm	ICP-MS 1000 mg/L Thulium Standard in %2 HNO3	%2 HNO3	250 ml	1000 mg/L
ICP500.TM1000.2NA	Tm	ICP-MS 1000 mg/L Thulium Standard in %2 HNO3	%2 HNO3	500 ml	1000 mg/L
ICP125.YB100.2NA	Yb	ICP-MS 100 mg/L Ytterbium Standard in %2 HNO3	%2 HNO3	125 ml	100 mg/L
ICP250.YB1000.2NA	Yb	ICP-MS 1000 mg/L Ytterbium Standard in %2 HNO3	%2 HNO3	250 ml	1000 mg/L
ICP500.YB1000.2NA	Yb	ICP-MS 1000 mg/L Ytterbium Standard in %2 HNO3	%2 HNO3	500 ml	1000 mg/L
ICP125.IR100.2CA	Ir	ICP-MS 100 mg/L Iridium Standard in %2 HCl	%2 HCl	125 ml	100 mg/L
ICP250.IR1000.2CA	Ir	ICP-MS 1000 mg/L Iridium Standard in %2 HCl	%2 HCl	250 ml	1000 mg/L
ICP500.IR1000.2CA	Ir	ICP-MS 1000 mg/L Iridium Standard in %2 HCl	%2 HCl	500 ml	1000 mg/L
ICP125.GA100.2NA	Ga	ICP-MS 100 mg/L Gallium Standard in %2 HNO3	%2 HNO3	125 ml	100 mg/L
ICP250.GA1000.2NA	Ga	ICP-MS 1000 mg/L Gallium Standard in %2 HNO3	%2 HNO3	250 ml	1000 mg/L
ICP500.GA1000.2NA	Ga	ICP-MS 1000 mg/L Gallium Standard in %2 HNO3	%2 HNO3	500 ml	1000 mg/L

ICP-MS Single Standards

ICP-MS Single Standards

ICP250.RE1000.2NA	Re	ICP-MS 1000 mg/L Rhenium Standard in %2 HNO3	%2 HNO3	250 ml	1000 mg/L
ICP500.RE1000.2NA	Re	ICP-MS 1000 mg/L Rhenium Standard in %2 HNO3	%2 HNO3	500 ml	1000 mg/L
ICP125.HF100.2NA05HF	Hf	ICP-MS 100 mg/L Hafnium Standard in %2 HNO3 + %0.5 HF	%2 HNO3 + Tr HF	125 ml	100 mg/L
ICP250.HF1000.2NA05HF	Hf	ICP-MS 1000 mg/L Hafnium Standard in %2 HNO3 + %0.5 HF	%2 HNO3 + Tr HF	250 ml	1000 mg/L
ICP500.HF1000.2NA05HF	Hf	ICP-MS 1000 mg/L Hafnium Standard in %2 HNO3 + %0.5 HF	%2 HNO3 + Tr HF	500 ml	1000 mg/L
ICP125.SN100.1NA05HF	Sn	ICP-MS 100 mg/L Tin Standard in %1 HNO3 + %0.5 HF	%1 HNO3 + Tr HF	125 ml	100 mg/L
ICP250.SN1000.1NA05HF	Sn	ICP-MS 1000 mg/L Tin Standard in %1 HNO3 + %0.5 HF	%1 HNO3 + Tr HF	250 ml	1000 mg/L
ICP500.SN1000.1NA05HF	Sn	ICP-MS 1000 mg/L Tin Standard in %1 HNO3 + %0.5 HF	%1 HNO3 + Tr HF	500 ml	1000 mg/L
ICP125.NB100.2NA05HF	Nb	ICP-MS 100 mg/L Niobium Standard in %2 HNO3 + %0.5 HF	%2 HNO3 + Tr HF	125 ml	100 mg/L
ICP250.NB1000.2NA05HF	Nb	ICP-MS 1000 mg/L Niobium Standard in %2 HNO3 + %0.5 HF	%2 HNO3 + Tr HF	250 ml	1000 mg/L
ICP500.NB1000.2NA05HF	Nb	ICP-MS 1000 mg/L Niobium Standard in %2 HNO3 + %0.5 HF	%2 HNO3 + Tr HF	500 ml	1000 mg/L
ICP125.TA100.2NA05HF	Ta	ICP-MS 100 mg/L Tantalum Standard in %2 HNO3 + %0.5 HF	%2 HNO3 + Tr HF	125 ml	100 mg/L
ICP250.TA1000.2NA05HF	Ta	ICP-MS 1000 mg/L Tantalum Standard in %2 HNO3 + %0.5 HF	%2 HNO3 + Tr HF	250 ml	1000 mg/L
ICP500.TA1000.2NA05HF	Ta	ICP-MS 1000 mg/L Tantalum Standard in %2 HNO3 + %0.5 HF	%2 HNO3 + Tr HF	500 ml	1000 mg/L
ICP125.ZR100.2NA05HF	Zr	ICP-MS 100 mg/L Zirconium Standard in %2 HNO3 + %0.5 HF	%2 HNO3 + Tr HF	125 ml	100 mg/L
ICP250.ZR1000.2NA05HF	Zr	ICP-MS 1000 mg/L Zirconium Standard in %2 HNO3 + %0.5 HF	%2 HNO3 + Tr HF	250 ml	1000 mg/L
ICP500.ZR1000.2NA05HF	Zr	ICP-MS 1000 mg/L Zirconium Standard in %2 HNO3 + %0.5 HF	%2 HNO3 + Tr HF	500 ml	1000 mg/L
ICP125.SB100.5NA1HF	Sb	ICP-MS 100 mg/L Antimony Standard in %5 HNO3 + %1 HF	%5 HNO3 + Tr HF	125 ml	100 mg/L
ICP250.SB1000.5NA1HF	Sb	ICP-MS 1000 mg/L Antimony Standard in %5 HNO3 + %1 HF	%5 HNO3 + Tr HF	250 ml	1000 mg/L
ICP500.SB1000.5NA1HF	Sb	ICP-MS 1000 mg/L Antimony Standard in %5 HNO3 + %1 HF	%5 HNO3 + Tr HF	500 ml	1000 mg/L
ICP125.TE100.2NA	Te	ICP-MS 100 mg/L Tellurium Standard in %2 HNO3	%2 HNO3	125 ml	100 mg/L
ICP250.TE1000.2NA	Te	ICP-MS 1000 mg/L Tellurium Standard in %2 HNO3	%2 HNO3	250 ml	1000 mg/L
ICP500.TE1000.2NA	Te	ICP-MS 1000 mg/L Tellurium Standard in %2 HNO3	%2 HNO3	500 ml	1000 mg/L
ICP125.TH100.2NA	Th	ICP-MS 100 mg/L Thorium Standard in %2 HNO3	%2 HNO3	125 ml	100 mg/L
ICP250.TH1000.2NA	Th	ICP-MS 1000 mg/L Thorium Standard in %2 HNO3	%2 HNO3	250 ml	1000 mg/L
ICP500.TH1000.2NA	Th	ICP-MS 1000 mg/L Thorium Standard in %2 HNO3	%2 HNO3	500 ml	1000 mg/L
ICP125.SM100.2NA	Sm	ICP-MS 100 mg/L Samarium Standard in %2 HNO3	%2 HNO3	125 ml	100 mg/L
ICP250.SM1000.2NA	Sm	ICP-MS 1000 mg/L Samarium Standard in %2 HNO3	%2 HNO3	250 ml	1000 mg/L
ICP500.SM1000.2NA	Sm	ICP-MS 1000 mg/L Samarium Standard in %2 HNO3	%2 HNO3	500 ml	1000 mg/L
ICP125.ND100.2NA	Nd	ICP-MS 100 mg/L Neodymium Standard in %2 HNO3	%2 HNO3	125 ml	100 mg/L
ICP250.ND1000.2NA	Nd	ICP-MS 1000 mg/L Neodymium Standard in %2 HNO3	%2 HNO3	250 ml	1000 mg/L
ICP500.ND1000.2NA	Nd	ICP-MS 1000 mg/L Neodymium Standard in %2 HNO3	%2 HNO3	500 ml	1000 mg/L
ICP125.MO100.W	Mo	ICP-MS 100 mg/L Molybdenum Standard in Water	Water	125 ml	100 mg/L
ICP250.MO1000.W	Mo	ICP-MS 1000 mg/L Molybdenum Standard in Water	Water	250 ml	1000 mg/L
ICP500.MO1000.W	Mo	ICP-MS 1000 mg/L Molybdenum Standard in Water	Water	500 ml	1000 mg/L
ICP125.P100.W	P	ICP-MS 100 mg/L Phosphorus Standard in Water	Water	125 ml	100 mg/L
ICP250.P1000.W	P	ICP-MS 1000 mg/L Phosphorus Standard in Water	Water	250 ml	1000 mg/L
ICP500.P1000.W	P	ICP-MS 1000 mg/L Phosphorus Standard in Water	Water	500 ml	1000 mg/L
ICP125.S100.W	S	ICP-MS 100 mg/L Sulphur Standard in Water	Water	125 ml	100 mg/L
ICP250.S1000.W	S	ICP-MS 1000 mg/L Sulphur Standard in Water	Water	250 ml	1000 mg/L

ICP-MS Single Standards

ICP-MS Single Standards

ICP500.S1000.W	S	ICP-MS 1000 mg/L Sulphur Standard in Water	Water	500 ml	1000 mg/L
ICP125.AG1000.5NA	Ag	ICP-MS 1000 mg/L Silver Standard in %5 HNO ₃	%5 HNO ₃	125 ml	1000 mg/L
ICP125.AL1000.1NA	Al	ICP-MS 1000 mg/L Aluminium Standard in %1 HNO ₃	%1 HNO ₃	125 ml	1000 mg/L
ICP125.AS1000.2CA	As	ICP-MS 1000 mg/L Arsenic Standard in %2 HCl	%2 HCl	125 ml	1000 mg/L
ICP125.AU1000.2CA	Au	ICP-MS 1000 mg/L Gold Standard in %2 HCl	%2 HCl	125 ml	1000 mg/L
ICP125.B1000.01NA	B	ICP-MS 1000 mg/L Boron Standard in %0.1HNO ₃	%0.1HNO ₃	125 ml	1000 mg/L
ICP125.BA1000.2NA	Ba	ICP-MS 1000 mg/L Barium Standard in %2 HNO ₃	%2 HNO ₃	125 ml	1000 mg/L
ICP125.BE1000.3NA	Be	ICP-MS 1000 mg/L Beryllium Standard in %3 HNO ₃	%3 HNO ₃	125 ml	1000 mg/L
ICP125.BI1000.3NA	Bi	ICP-MS 1000 mg/L Bismuth Standard in %3 HNO ₃	%3 HNO ₃	125 ml	1000 mg/L
ICP125.CA1000.3NA	Ca	ICP-MS 1000 mg/L Calcium Standard in %0.1 HNO ₃	%0.1 HNO ₃	125 ml	1000 mg/L
ICP125.CD1000.3NA	Cd	ICP-MS 1000 mg/L Cadmium Standard in %3 HNO ₃	%3 HNO ₃	125 ml	1000 mg/L
ICP125.CE1000.3NA	Ce	ICP-MS 1000 mg/L Cerium Standard in %3 HNO ₃	%3 HNO ₃	125 ml	1000 mg/L
ICP125.CO1000.2NA	Co	ICP-MS 1000 mg/L Cobalt Standard in %2 HNO ₃	%2 HNO ₃	125 ml	1000 mg/L
ICP125.CS1000.2NA	Cs	ICP-MS 1000 mg/L Cesium Standard in %2 HNO ₃	%2 HNO ₃	125 ml	1000 mg/L
ICP125.CU1000.2NA	Cu	ICP-MS 1000 mg/L Copper Standard in %2 HNO ₃	%2 HNO ₃	125 ml	1000 mg/L
ICP125.DY1000.2NA	Dy	ICP-MS 1000 mg/L Dysprosium Standard in %2 HNO ₃	%2 HNO ₃	125 ml	1000 mg/L
ICP125.EU1000.2NA	Eu	ICP-MS 1000 mg/L Europium Standard in %2 HNO ₃	%2 HNO ₃	125 ml	1000 mg/L
ICP125.FE1000.3NA	Fe	ICP-MS 1000 mg/L Iron Standard in %3 HNO ₃	%3 HNO ₃	125 ml	1000 mg/L
ICP125.GA1000.2NA	Ga	ICP-MS 1000 mg/L Gallium Standard in %2 HNO ₃	%2 HNO ₃	125 ml	1000 mg/L
ICP125.GD1000.2NA	Gd	ICP-MS 1000 mg/L Gadolinium Standard in %2 HNO ₃	%2 HNO ₃	125 ml	1000 mg/L
ICP125.GE10.1NA	Ge	ICP-MS 10 mg/L Germanium Standard in %1 HNO ₃	%1 HNO ₃	125 ml	10 mg/L
ICP125.HF1000.2NA01HF	Hf	ICP-MS 1000 mg/L Hafnium Standard in %2 HNO ₃ + % 0.1 HF	%2 HNO ₃ + Tr HF	125 ml	1000 mg/L
ICP125.HG1000.2CA	Hg	ICP-MS 1000 mg/L Mercury Standard in %2 HCl	%2 HCl	125 ml	1000 mg/L
ICP125.HO1000.2NA	Ho	ICP-MS 1000 mg/L Holmium Standard in %2 HNO ₃	%2 HNO ₃	125 ml	1000 mg/L
ICP125.IN1000.3NA	In	ICP-MS 1000 mg/L Indium Standard in %3 HNO ₃	%3 HNO ₃	125 ml	1000 mg/L
ICP125.IR1000.2CA	Ir	ICP-MS 1000 mg/L Iridium Standard in %2 HCl	%2 HCl	125 ml	1000 mg/L
ICP125.K1000.01NA	K	ICP-MS 1000 mg/L Potassium Standard in %0.1 HNO ₃	%0.1 HNO ₃	125 ml	1000 mg/L
ICP125.LA1000.2NA	La	ICP-MS 1000 mg/L Lanthanum Standard in %2 HNO ₃	%2 HNO ₃	125 ml	1000 mg/L
ICP125.LI1000.01NA	Li	ICP-MS 1000 mg/L Lithium Standard in %0.1 HNO ₃	%0.1 HNO ₃	125 ml	1000 mg/L
ICP125.LU1000.3NA	Lu	ICP-MS 1000 mg/L Lutetium Standard in %3 HNO ₃	%3 HNO ₃	125 ml	1000 mg/L
ICP125.MG1000.01NA	Mg	ICP-MS 1000 mg/L Magnesium Standard in %0.1 HNO ₃	%0.1 HNO ₃	125 ml	1000 mg/L
ICP125.MN1000.3NA	Mn	ICP-MS 1000 mg/L Manganese Standard in %3 HNO ₃	%3 HNO ₃	125 ml	1000 mg/L
ICP125.MO1000.W	Mo	ICP-MS 1000 mg/L Molybdenum Standard in Water	Water	125 ml	1000 mg/L
ICP125.NA1000.01NA	Na	ICP-MS 1000 mg/L Sodium Standard in %0.1 HNO ₃	%0.1 HNO ₃	125 ml	1000 mg/L
ICP125.NB1000.2NA05HF	Nb	ICP-MS 1000 mg/L Niobium Standard in %2 HNO ₃ + %0.5 HF	%2 HNO ₃ + Tr HF	125 ml	1000 mg/L
ICP125.ND1000.2NA	Nd	ICP-MS 1000 mg/L Neodymium Standard in %2 HNO ₃	%2 HNO ₃	125 ml	1000 mg/L
ICP125.NH41000.01NA	NH ₄	ICP-MS 1000 mg/L Ammonium Standard in %0.1 HNO ₃	%0.1 HNO ₃	125 ml	1000 mg/L
ICP125.NI1000.2NA	Ni	ICP-MS 1000 mg/L Nickel Standard in %2 HNO ₃	%2 HNO ₃	125 ml	1000 mg/L
ICP125.P1000.W	P	ICP-MS 1000 mg/L Phosphorus Standard in Water	Water	125 ml	1000 mg/L
ICP125.PB1000.05NA	Pb	ICP-MS 1000 mg/L Lead Standard in %0.5 HNO ₃	%0.5 HNO ₃	125 ml	1000 mg/L

ICP-MS Single Standards

ICP-MS Single Standards

ICP125.PD1000.05NA	Pd	ICP-MS 1000 mg/L Palladium Standard in %0.5 HNO ₃	%0.5 HNO ₃	125 ml	1000 mg/L
ICP125.PR1000.2NA	Pr	ICP-MS 1000 mg/L Praseodymium Standard in %2 HNO ₃	%2 HNO ₃	125 ml	1000 mg/L
ICP125.PT1000.2CA	Pt	ICP-MS 1000 mg/L Platinum Standard in %2 HCl	%2 HCl	125 ml	1000 mg/L
ICP125.RB1000.2NA	Rb	ICP-MS 1000 mg/L Rubidium Standard in %2 HNO ₃	%2 HNO ₃	125 ml	1000 mg/L
ICP125.RE1000.2NA	Re	ICP-MS 1000 mg/L Rhenium Standard in %2 HNO ₃	%2 HNO ₃	125 ml	1000 mg/L
ICP125.RH10.3NA	Rh	ICP-MS 10 mg/L Rhodium Standard in %3 HNO ₃	%3 HNO ₃	125 ml	10 mg/L
ICP125.RU1000.2CA	Ru	ICP-MS 1000 mg/L Ruthenium Standard in %2 HCl	%2 HCl	125 ml	1000 mg/L
ICP125.S1000.W	S	ICP-MS 1000 mg/L Sulphur Standard in Water	Water	125 ml	1000 mg/L
ICP125.SB1000.5NA1HF	Sb	ICP-MS 1000 mg/L Antimony Standard in %5 HNO ₃ + %1 HF	%5 HNO ₃ + Tr HF	125 ml	1000 mg/L
ICP125.SC10.2NA	Sc	ICP-MS 10 mg/L Scandium Standard in %2 HNO ₃	%2 HNO ₃	125 ml	10 mg/L
ICP125.SE1000.2NA	Se	ICP-MS 1000 mg/L Selenium Standard in %2 HNO ₃	%2 HNO ₃	125 ml	1000 mg/L
ICP125.SI1000.3NA	Si	ICP-MS 1000 mg/L Silicon Standard in %3 HNO ₃	%3 HNO ₃	125 ml	1000 mg/L
ICP125.SM1000.2NA	Sm	ICP-MS 1000 mg/L Samarium Standard in %2 HNO ₃	%2 HNO ₃	125 ml	1000 mg/L
ICP125.SN1000.1NA05HF	Sn	ICP-MS 1000 mg/L Tin Standard in %1 HNO ₃ + %0.5 HF	%1 HNO ₃ + Tr HF	125 ml	1000 mg/L
ICP125.SR1000.2NA	Sr	ICP-MS 1000 mg/L Strontium Standard in %2 HNO ₃	%2 HNO ₃	125 ml	1000 mg/L
ICP125.TA1000.2NA05HF	Ta	ICP-MS 1000 mg/L Tantalum Standard in %2 HNO ₃ + %0.5 HF	%2 HNO ₃ + Tr HF	125 ml	1000 mg/L
ICP125.TB1000.3NA	Tb	ICP-MS 1000 mg/L Terbium Standard in %3 HNO ₃	%3 HNO ₃	125 ml	1000 mg/L
ICP125.TE1000.2NA	Te	ICP-MS 1000 mg/L Tellurium Standard in %2 HNO ₃	%2 HNO ₃	125 ml	1000 mg/L
ICP125.TH1000.2NA	Th	ICP-MS 1000 mg/L Thorium Standard in %2 HNO ₃	%2 HNO ₃	125 ml	1000 mg/L
ICP125.TI1000.3NA	Ti	ICP-MS 1000 mg/L Titanium Standard in %3 HNO ₃	%3 HNO ₃	125 ml	1000 mg/L
ICP125.TI1000.2NA	Tl	ICP-MS 1000 mg/L Thallium Standard in %2 HNO ₃	%2 HNO ₃	125 ml	1000 mg/L
ICP125.TM1000.2NA	Tm	ICP-MS 1000 mg/L Thulium Standard in %2 HNO ₃	%2 HNO ₃	125 ml	1000 mg/L
ICP125.U1000.2NA	U	ICP-MS 1000 mg/L Uranium Standard in %2 HNO ₃	%2 HNO ₃	125 ml	1000 mg/L
ICP125.V1000.2NA	V	ICP-MS 1000 mg/L Vanadium Standard in %2 HNO ₃	%2 HNO ₃	125 ml	1000 mg/L
ICP125.Y1000.2NA	Y	ICP-MS 1000 mg/L Yttrium Standard in %2 HNO ₃	%2 HNO ₃	125 ml	1000 mg/L
ICP125.YB1000.2NA	Yb	ICP-MS 1000 mg/L Ytterbium Standard in %2 HNO ₃	%2 HNO ₃	125 ml	1000 mg/L
ICP125.ZN1000.2NA	Zn	ICP-MS 1000 mg/L Zinc Standard in %2 HNO ₃	%2 HNO ₃	125 ml	1000 mg/L
ICP125.ZR1000.2NA05HF	Zr	ICP-MS 1000 mg/L Zirconium Standard in %2 HNO ₃ + %0.5 HF	%2 HNO ₃ + Tr HF	125 ml	1000 mg/L
ICP125.CR1000.2NA	Cr	ICP-MS 1000 mg/L Chromium Standard in %2 HNO ₃	%2 HNO ₃	125 ml	1000 mg/L
ICP125.ER1000.2NA	Er	ICP-MS 1000 mg/L Erbium Standard in %2 HNO ₃	%2 HNO ₃	125 ml	1000 mg/L
ICP125.CR31000.5NA	Cr	ICP-MS 1000 mg/L Chromium (III) Standard in %5 HNO ₃	%5 HNO ₃	125 ml	1000 mg/L
ICP250.CR31000.5NA	Cr	ICP-MS 1000 mg/L Chromium (III) Standard in %5 HNO ₃	%5 HNO ₃	250 ml	1000 mg/L
ICP500.CR31000.5N	Cr	ICP-MS 1000 mg/L Chromium (III) Standard in %5 HNO ₃	%5 HNO ₃	500 ml	1000 mg/L

We would like to offer tailor-made products for your special studies. Please contact us for custom-made reference standard requests.

ICP-MS Multi Element Standards

ICP-MS Multi Element Standards

Product Code:	ICP125.M28001	Matrix:	3% HNO ₃
Chapter:	125ml-Tuning Solution - 28001	Volume:	125 ml
Contents			
Be: 10, Co: 10, In: 10, Pb: 10, Mg: 10 (mg/l)			

Product Code:	ICP125.M28002	Matrix:	5% HNO ₃
Chapter:	125ml-Internal Standard - 28002	Volume:	125 ml
Contents			
Bi: 10, Ge: 10, In: 10, 6Li: 10, Sc: 10, Tb: 10, Y: 10 (mg/l)			

Product Code:	ICP125.M28003	Matrix:	2% HNO ₃
Chapter:	125ml-Calibration Standard - 28003	Volume:	125 ml
Contents			
10 mg/l [Al] Aluminium, 10 mg/l [As] Arsenic, 10 mg/l [Ba] Barium, 10 mg/l [Be] Beryllium, 10 mg/l [Cd] Cadmium, 10 mg/l [Cr] Chromium, 10 mg/l [Co] Cobalt, 10 mg/l [Cu] Copper, 10 mg/l [Pb] Lead, 10 mg/l [Mn] Manganese, 10 mg/l [Ni] Nickel, 10 mg/l [Se] Selenium, 10 mg/l [Ag] Silver, 10 mg/l [Tl] Thallium, 10 mg/l [Th] Thorium, 10 mg/l [U] Uranium, 10 mg/l [V] Vanadium, 10 mg/l [Zn] Zinc			

Product Code:	ICP125.M28004	Matrix:	5% HCl
Chapter:	125ml-Calibration Standard - 28004	Volume:	125 ml
Contents			
20 mg/l [Al] Aluminium, 20 mg/l [As] Arsenic, 2 mg/l [Ba] Barium, 1 mg/l [Be] Beryllium, 2 mg/l [B] Boron, 2 mg/l [Cd] Cadmium, 10 mg/l [Ca] Calcium, 2 mg/l [Cr] Chromium, 2 mg/l [Cu] Copper, 2 mg/l [Fe] Iron, 20 mg/l [Pb] Lead, 2 mg/l [Li] Lithium, 1 mg/l [Mg] Magnesium, 1 mg/l [Mn] Manganese, 5 mg/l [Hg] Mercury, 5 mg/l [Ni] Nickel, 10 mg/l [P] Phosphorus, 100 mg/l [K] Potassium, 1 mg/l [Sc] Scandium, 20 mg/l [Se] Selenium, 20 mg/l [Na] Sodium, 1 mg/l [Sr] Strontium, 20 mg/l [Te] Tellurium, 2 mg/l [Ti] Titanium, 1 mg/l [Y] Yttrium, 2 mg/l [Zn] Zinc			

ICP-MS Multi Element Standards

ICP-MS Multi Element Standards

Product Code:	ICP125.M28005	Matrix:	5% HCl
Chapter:	125ml-Precious Metals Standard - 28005	Volume:	125 ml
Contents			
10 mg/l [Au] Gold, 10 mg/l [Ir] Iridium, 10 mg/l [Os] Osmium, 10 mg/l [Pd] Palladium, 10 mg/l [Pt] Platinum, 10 mg/l [Re] Rhenium, 10 mg/l [Rh] Rhodium, 10 mg/l [Ru] Ruthenium			

Product Code:	ICP125.M28006	Matrix:	5% HNO ₃
Chapter:	125ml-Rare Earth Metal Standard - 28006	Volume:	125 ml
Contents			
10 ug/ml [Ce] Cerium, 10 ug/ml [Dy] Dysprosium, 10 ug/ml [Er] Erbium, 10 ug/ml [Eu] Europium, 10 ug/ml [Gd] Gadolinium, 10 ug/ml [Ho] Holmium, 10 ug/ml [La] Lanthanum, 10 ug/ml [Lu] Lutetium, 10 ug/ml [Nd] Neodymium, 10 ug/ml [Pr] Praseodymium, 10 ug/ml [Sm] Samarium, 10 ug/ml [Sc] Scandium, 10 ug/ml [Tb] Terbium, 10 ug/ml [Th] Thorium, 10 ug/ml [Tm] Thulium, 10 ug/ml [U] Uranium, 10 ug/ml [Yb] Ytterbium, 10 ug/ml [Y] Yttrium			

Product Code:	ICP125.M28007	Matrix:	5% HNO ₃
Chapter:	125ml-Calibration Standard - 28007	Volume:	125 ml
Contents			
100 mg/l [Be] Beryllium, 1000 mg/l [Fe] Iron, 1000 mg/l [Pb] Lead, 1000 mg/l [Mg] Magnesium, 200 mg/l [Ni] Nickel, 500 mg/l [Tl] Thallium			

Product Code:	ICP125.M28008	Matrix:	5% HNO ₃
Chapter:	125ml-Calibration Standard - 28008	Volume:	125 ml
Contents			
[Sb] Antimony 500 mg/l, [As] Arsenic 1000 mg/l, [Ba] Barium 100 mg/l, [B] Boron 100 mg/l, [Cd] Cadmium 200 mg/l, [Ca] Calcium 1000 mg/l, [Cu] Copper 200 mg/l, [Mn] Manganese 200 mg/l, [Se] Selenium 500 mg/l, [Ag] Silver 50 mg/l			

ICP-MS Multi Element Standards

ICP-MS Multi Element Standards

Product Code:	ICP125.M28009	Matrix:	5% HNO ₃
Chapter:	125ml-Calibration Standard - 28009	Volume:	125 ml
Contents			
[Al] Aluminium 1000 mg/l, [Cr] Chromium 500 mg/l, [Hg] Mercury 200 mg/l, [Zn] Zinc 500 mg/l			

Product Code:	ICP125.M28010	Matrix:	2% HNO ₃
Chapter:	125ml-Tuning Solution - 28010	Volume:	125 ml
Contents			
10 mg/l [Ba] Barium, 10 mg/l [Be] Beryllium, 10 mg/l [Ce] Cerium, 10 mg/l [Co] Cobalt, 10 mg/l [In] Indium, 10 mg/l [Pb] Lead, 10 mg/l [Mg] Magnesium, 10 mg/l [Tl] Thallium, 10 mg/l [Th] Thorium			

Product Code:	ICP125.M28011	Matrix:	5% HNO ₃
Chapter:	125ml-Internal Standard - 28011	Volume:	125 ml
Contents			
100 mg/l [Bi] Bismuth, 100 mg/l [In] Indium, 100 mg/l [6Li] Lithium isotope 6, 100 mg/l [Sc] Scandium, 100 mg/l [Tb] Terbium, 100 mg/l [Y] Yttrium			

Product Code:	ICP125.M28012	Matrix:	5% HNO ₃ / Tr HF
Chapter:	125ml-Refractory Elements Standard - 28012	Volume:	125 ml
Contents			
10 mg/l [Sb] Antimony, 10 mg/l [Ge] Germanium, 10 mg/l [Hf] Hafnium, 10 mg/l [Mo] Molybdenum, 10 mg/l [Nb] Niobium, 10 mg/l [Si] Silicon, 10 mg/l [Ta] Tantalum, 10 mg/l [Te] Tellurium, 10 mg/l [Sn] Tin, 10 mg/l [Ti] Titanium, 10 mg/l [W] Tungsten, 10 mg/l [Zr] Zirconium			

ICP-MS Multi Element Standards

ICP-MS Multi Element Standards

Product Code:	ICP125.M28013	Matrix:	10% HNO ₃
Chapter:	125ml-Transition Metals Standard - 28013	Volume:	125 ml
Contents			
100 mg/l [Cd] Cadmium, 100 mg/l [Cr] Chromium, 100 mg/l [Co] Cobalt, 100 mg/l [Cu] Copper, 100 mg/l [Fe] Iron, 100 mg/l [Pb] Lead, 100 mg/l [Mn] Manganese, 100 mg/l [Hg] Mercury, 100 mg/l [Ni] Nickel, 100 mg/l [Ag] Silver, 100 mg/l [Tl] Thallium, 100 mg/l [V] Vanadium, 100 mg/l [Zn] Zinc			

Product Code:	ICP125.M28014	Matrix:	5% HNO ₃ / Tr HF
Chapter:	125ml-Instrument Fortifying Standard - 28014	Volume:	125 ml
Contents			
[Sb] Antimony 20 mg/l, [Mo] Molybdenum 20 mg/l, [Si] Silicon 20 mg/l, [Sn] Tin 20 mg/l, [Ti] Titanium 20 mg/l			

Product Code:	ICP125.M28016	Matrix:	%5 HNO ₃
Chapter:	125ml-Calibration Standard - 28016	Volume:	125 ml
Contents			
100 mg/l [Al] Aluminium, 100 mg/l [As] Arsenic, 100 mg/l [Ba] Barium, 100 mg/l [Be] Beryllium, 100 mg/l [Bi] Bismuth, 100 mg/l [B] Boron, 100 mg/l [Cd] Cadmium, 100 mg/l [Ca] Calcium, 100 mg/l [Cr] Chromium, 100 mg/l [Co] Cobalt, 100 mg/l [Cu] Copper, 100 mg/l [Fe] Iron, 100 mg/l [Pb] Lead, 100 mg/l [Li] Lithium, 100 mg/l [Mg] Magnesium, 100 mg/l [Mn] Manganese, 100 mg/l [Mo] Molybdenum, 100 mg/l [Ni] Nickel, 100 mg/l [K] Potassium, 100 mg/l [Se] Selenium, 100 mg/l [Si] Silicon, 100 mg/l [Na] Sodium, 100 mg/l [Sr] Strontium, 100 mg/l [Tl] Thallium, 100 mg/l [Ti] Titanium, 100 mg/l [V] Vanadium, 100 mg/l [Zn] Zinc			



We manufacture and distribute our products for laboratories around the world. We arrange the most suitable shipment terms for your orders.

ICP-MS Multi Element Standards

ICP-MS Multi Element Standards

Product Code:	ICP125.M28017	Matrix:	%5 HNO ₃ + tr HF
Chapter:	125ml-Calibration Standard - 28017	Volume:	125 ml
Contents			
100 mg/l [P] Phosphorus, 5 mg/l [Ag] Silver, 20 mg/l [Sb] Antimony, 20 mg/l [As] Arsenic, 20 mg/l [Ba] Barium, 20 mg/l [Be] Beryllium, 20 mg/l [Cd] Cadmium, 20 mg/l [Cr] Chromium, 20 mg/l [Co] Cobalt, 20 mg/l [Cu] Copper, 20 mg/l [Pb] Lead, 20 mg/l [Mn] Manganese, 20 mg/l [Ni] Nickel, 20 mg/l [Se] Selenium, 20 mg/l [Zn] Zinc, 20 mg/l [V] Vanadium, 20 mg/l [Tl] Thallium, 20 mg/l [Fe] Iron			

Product Code:	ICP250.M28001	Matrix:	3% HNO ₃
Chapter:	250ml-Tuning Solution - 28001	Volume:	250 ml
Contents			
Be: 10, Co: 10, In: 10, Pb: 10, Mg: 10 (mg/l)			

Product Code:	ICP250.M28002	Matrix:	5% HNO ₃
Chapter:	250ml-Internal Standard - 28002	Volume:	250 ml
Contents			
Bi: 10, Ge: 10, In: 10, 6Li: 10, Sc: 10, Tb: 10, Y: 10 (mg/l)			

ICP-MS Multi Element Standards

ICP-MS Multi Element Standards

Product Code:	ICP250.M28003	Matrix:	2% HNO ₃
Chapter:	250ml-Calibration Standard - 28003	Volume:	250 ml
Contents			
10 mg/l [Al] Aluminium, 10 mg/l [As] Arsenic, 10 mg/l [Ba] Barium, 10 mg/l [Be] Beryllium, 10 mg/l [Cd] Cadmium, 10 mg/l [Cr] Chromium, 10 mg/l [Co] Cobalt, 10 mg/l [Cu] Copper, 10 mg/l [Pb] Lead, 10 mg/l [Mn] Manganese, 10 mg/l [Ni] Nickel, 10 mg/l [Se] Selenium, 10 mg/l [Ag] Silver, 10 mg/l [Tl] Thallium, 10 mg/l [Th] Thorium, 10 mg/l [U] Uranium, 10 mg/l [V] Vanadium, 10 mg/l [Zn] Zinc			

Product Code:	ICP250.M28004	Matrix:	5% HCl
Chapter:	250ml-Calibration Standard - 28004	Volume:	250 ml
Contents			
20 mg/l [Al] Aluminium, 20 mg/l [As] Arsenic, 2 mg/l [Ba] Barium, 1 mg/l [Be] Beryllium, 2 mg/l [B] Boron, 2 mg/l [Cd] Cadmium, 10 mg/l [Ca] Calcium, 2 mg/l [Cr] Chromium, 2 mg/l [Cu] Copper, 2 mg/l [Fe] Iron, 20 mg/l [Pb] Lead, 2 mg/l [Li] Lithium, 1 mg/l [Mg] Magnesium, 1 mg/l [Mn] Manganese, 5 mg/l [Hg] Mercury, 5 mg/l [Ni] Nickel, 10 mg/l [P] Phosphorus, 100 mg/l [K] Potassium, 1 mg/l [Sc] Scandium, 20 mg/l [Se] Selenium, 20 mg/l [Na] Sodium, 1 mg/l [Sr] Strontium, 20 mg/l [Te] Tellurium, 2 mg/l [Ti] Titanium, 1 mg/l [Y] Yttrium, 2 mg/l [Zn] Zinc			

Product Code:	ICP250.M28005	Matrix:	5% HCl
Chapter:	250ml-Precious Metals Standard - 28005	Volume:	250 ml
Contents			
10 mg/l [Au] Gold, 10 mg/l [Ir] Iridium, 10 mg/l [Os] Osmium, 10 mg/l [Pd] Palladium, 10 mg/l [Pt] Platinum, 10 mg/l [Re] Rhenium, 10 mg/l [Rh] Rhodium, 10 mg/l [Ru] Ruthenium			

Product Code:	ICP250.M28006	Matrix:	5% HNO ₃
Chapter:	250ml-Rare Earth Metal Standard - 28006	Volume:	250 ml
Contents			
10 ug/ml [Ce] Cerium, 10 ug/ml [Dy] Dysprosium, 10 ug/ml [Er] Erbium, 10 ug/ml [Eu] Europium, 10 ug/ml [Gd] Gadolinium, 10 ug/ml [Ho] Holmium, 10 ug/ml [La] Lanthanum, 10 ug/ml [Lu] Lutetium, 10 ug/ml [Nd] Neodymium, 10 ug/ml [Pr] Praseodymium, 10 ug/ml [Sm] Samarium, 10 ug/ml [Sc] Scandium, 10 ug/ml [Tb] Terbium, 10 ug/ml [Th] Thorium, 10 ug/ml [Tm] Thulium, 10 ug/ml [U] Uranium, 10 ug/ml [Yb] Ytterbium, 10 ug/ml [Y] Yttrium			

ICP-MS Multi Element Standards

ICP-MS Multi Element Standards

Product Code:	ICP250.M28007	Matrix:	5% HNO ₃
Chapter:	250ml-Calibration Standard - 28007	Volume:	250 ml
Contents			
100 mg/l [Be] Beryllium, 1000 mg/l [Fe] Iron, 1000 mg/l [Pb] Lead, 1000 mg/l [Mg] Magnesium, 200 mg/l [Ni] Nickel, 500 mg/l [Tl] Thallium			

Product Code:	ICP250.M28008	Matrix:	5% HNO ₃
Chapter:	250ml-Calibration Standard - 28008	Volume:	250 ml
Contents			
[Sb] Antimony 500 mg/l, [As] Arsenic 1000 mg/l, [Ba] Barium 100 mg/l, [B] Boron 100 mg/l, [Cd] Cadmium 200 mg/l, [Ca] Calcium 1000 mg/l, [Cu] Copper 200 mg/l, [Mn] Manganese 200 mg/l, [Se] Selenium 500 mg/l, [Ag] Silver 50 mg/l			

Product Code:	ICP250.M28009	Matrix:	5% HNO ₃
Chapter:	250ml-Calibration Standard - 28009	Volume:	250 ml
Contents			
[Al] Aluminium 1000 mg/l, [Cr] Chromium 500 mg/l, [Hg] Mercury 200 mg/l, [Zn] Zinc 500 mg/l			

Product Code:	ICP250.M28010	Matrix:	2% HNO ₃
Chapter:	250ml-Tuning Solution - 28010	Volume:	250 ml
Contents			
10 mg/l [Ba] Barium, 10 mg/l [Be] Beryllium, 10 mg/l [Ce] Cerium, 10 mg/l [Co] Cobalt, 10 mg/l [In] Indium, 10 mg/l [Pb] Lead, 10 mg/l [Mg] Magnesium, 10 mg/l [Tl] Thallium, 10 mg/l [Th] Thorium			

ICP-MS Multi Element Standards

ICP-MS Multi Element Standards

Product Code:	ICP250.M28011	Matrix:	5% HNO ₃
Chapter:	250ml-Internal Standard - 28011	Volume:	250 ml
Contents			
100 mg/l [Bi] Bismuth, 100 mg/l [In] Indium, 100 mg/l [6Li] Lithium isotope 6, 100 mg/l [Sc] Scandium, 100 mg/l [Tb] Terbium, 100 mg/l [Y] Yttrium			

Product Code:	ICP250.M28012	Matrix:	5% HNO ₃ / Tr HF
Chapter:	250ml-Refractory Elements Standard - 28012	Volume:	250 ml
Contents			
10 mg/l [Sb] Antimony, 10 mg/l [Ge] Germanium, 10 mg/l [Hf] Hafnium, 10 mg/l [Mo] Molybdenum, 10 mg/l [Nb] Niobium, 10 mg/l [Si] Silicon, 10 mg/l [Ta] Tantalum, 10 mg/l [Te] Tellurium, 10 mg/l [Sn] Tin, 10 mg/l [Ti] Titanium, 10 mg/l [W] Tungsten, 10 mg/l [Zr] Zirconium			

Product Code:	ICP250.M28013	Matrix:	10% HNO ₃
Chapter:	250ml-Transition Metals Standard - 28013	Volume:	250 ml
Contents			
100 mg/l [Cd] Cadmium, 100 mg/l [Cr] Chromium, 100 mg/l [Co] Cobalt, 100 mg/l [Cu] Copper, 100 mg/l [Fe] Iron, 100 mg/l [Pb] Lead, 100 mg/l [Mn] Manganese, 100 mg/l [Hg] Mercury, 100 mg/l [Ni] Nickel, 100 mg/l [Ag] Silver, 100 mg/l [Tl] Thallium, 100 mg/l [V] Vanadium, 100 mg/l [Zn] Zinc			

Product Code:	ICP250.M28014	Matrix:	5% HNO ₃ / Tr HF
Chapter:	250ml-Instrument Fortifying Standard - 28014	Volume:	250 ml
Contents			
[Sb] Antimony 20 mg/l, [Mo] Molybdenum 20 mg/l, [Si] Silicon 20 mg/l, [Sn] Tin 20 mg/l, [Ti] Titanium 20 mg/l			


Product Code:	ICP250.M28016	Matrix:	%5 HNO ₃
Chapter:	250ml-Calibration Standard - 28016	Volume:	250 ml
Contents			
100 mg/l [Al] Aluminium, 100 mg/l [As] Arsenic, 100 mg/l [Ba] Barium, 100 mg/l [Be] Beryllium, 100 mg/l [Bi] Bismuth, 100 mg/l [B] Boron, 100 mg/l [Cd] Cadmium, 100 mg/l [Ca] Calcium, 100 mg/l [Cr] Chromium, 100 mg/l [Co] Cobalt, 100 mg/l [Cu] Copper, 100 mg/l [Fe] Iron, 100 mg/l [Pb] Lead, 100 mg/l [Li] Lithium, 100 mg/l [Mg] Magnesium, 100 mg/l [Mn] Manganese, 100 mg/l [Mo] Molybdenum, 100 mg/l [Ni] Nickel, 100 mg/l [K] Potassium, 100 mg/l [Se] Selenium, 100 mg/l [Si] Silicon, 100 mg/l [Na] Sodium, 100 mg/l [Sr] Strontium, 100 mg/l [Tl] Thallium, 100 mg/l [Ti] Titanium, 100 mg/l [V] Vanadium, 100 mg/l [Zn] Zinc			

Product Code:	ICP250.M28017	Matrix:	%5 HNO ₃ + tr HF
Chapter:	250ml-Calibration Standard - 28017	Volume:	250 ml
Contents			
100 mg/l [P] Phosphorus, 5 mg/l [Ag] Silver, 20 mg/l [Sb] Antimony, 20 mg/l [As] Arsenic, 20 mg/l [Ba] Barium, 20 mg/l [Be] Beryllium, 20 mg/l [Cd] Cadmium, 20 mg/l [Cr] Chromium, 20 mg/l [Co] Cobalt, 20 mg/l [Cu] Copper, 20 mg/l [Pb] Lead, 20 mg/l [Mn] Manganese, 20 mg/l [Ni] Nickel, 20 mg/l [Se] Selenium, 20 mg/l [Zn] Zinc, 20 mg/l [V] Vanadium, 20 mg/l [Tl] Thallium, 20 mg/l [Fe] Iron			

ICP-MS Blank Solutions

Product Code	Symbol	Description	Matrix	Volume	Value
BK500.HNO3.5		500 ml Sample Preparation Blank Solution- HNO ₃ 5%	Water	500 ml	%5 HNO ₃
BK1000.HNO3.5		1000 ml Sample Preparation Blank Solution - HNO ₃ 5%	Water	1000 ml	%5 HNO ₃
BK500.H2O.100		500 ml Blank Solution - High Purity Water for Chromatography	Water	500 ml	xx
BK1000.H2O.100		1000 ml Blank Solution - High Purity Water	Water	1000 ml	xx

ION CHROMATOGRAPHY STANDARDS



Ion Chromatography Single Ion Standards

Ion Chromatography Single Ion Standards

Product Code	Symbol	Description	Matrix	Volume	Value
IC125.SO41000.W	SO4	1000 mg/L-125ml-Sulfate	Water	125 ml	1000 mg/L
IC250.SO41000.W	SO4	1000 mg/L-250ml-Sulfate	Water	250 ml	1000 mg/L
IC500.SO41000.W	SO4	1000 mg/L-500ml-Sulfate	Water	500 ml	1000 mg/L
IC125.K1000.W	K	1000 mg/L-125ml-Potassium	Water	125 ml	1000 mg/L
IC250.K1000.W	K	1000 mg/L-250ml-Potassium	Water	250 ml	1000 mg/L
IC500.K1000.W	K	1000 mg/L-500ml-Potassium	Water	500 ml	1000 mg/L
IC125.NO31000.W	NO3	1000 mg/L-125ml-Nitrate	Water	125 ml	1000 mg/L
IC250.NO31000.W	NO3	1000 mg/L-250ml-Nitrate	Water	250 ml	1000 mg/L
IC500.NO31000.W	NO3	1000 mg/L-500ml-Nitrate	Water	500 ml	1000 mg/L
IC125.NO21000.W	NO2	1000 mg/L-125ml-Nitrite	Water	125 ml	1000 mg/L
IC250.NO21000.W	NO2	1000 mg/L-250ml-Nitrite	Water	250 ml	1000 mg/L
IC500.NO21000.W	NO2	1000 mg/L-500ml-Nitrite	Water	500 ml	1000 mg/L
IC125.F1000.W	F	1000 mg/L-125ml-Fluoride	Water	125 ml	1000 mg/L
IC250.F1000.W	F	1000 mg/L-250ml-Fluoride	Water	250 ml	1000 mg/L
IC500.F1000.W	F	1000 mg/L-500ml-Fluoride	Water	500 ml	1000 mg/L
IC125.BRO31000.W	BrO3	1000 mg/L-125ml-Bromate	Water	125 ml	1000 mg/L
IC250.BRO31000.W	BrO3	1000 mg/L-250ml-Bromate	Water	250 ml	1000 mg/L
IC500.BRO31000.W	BrO3	1000 mg/L-500ml-Bromate	Water	500 ml	1000 mg/L
IC125.BR1000.W	Br	1000 mg/L-125ml-Bromide	Water	125 ml	1000 mg/L
IC250.BR1000.W	Br	1000 mg/L-250ml-Bromide	Water	250 ml	1000 mg/L
IC500.BR1000.W	Br	1000 mg/L-500ml-Bromide	Water	500 ml	1000 mg/L
IC125.PO41000.W	PO4	1000 mg/L-125ml-Phosphate	Water	125 ml	1000 mg/L
IC250.PO41000.W	PO4	1000 mg/L-250ml-Phosphate	Water	250 ml	1000 mg/L
IC500.PO41000.W	PO4	1000 mg/L-500ml-Phosphate	Water	500 ml	1000 mg/L
IC125.NH41000.W	NH4	1000 mg/L-125ml-Ammonium	Water	125 ml	1000 mg/L
IC250.NH41000.W	NH4	1000 mg/L-250ml-Ammonium	Water	250 ml	1000 mg/L
IC500.NH41000.W	NH4	1000 mg/L-500ml-Ammonium	Water	500 ml	1000 mg/L
IC125.Li1000.W	Li	1000 mg/L -125ml- Lithium	Water	125 ml	1000 mg/L
IC250.Li1000.W	Li	1000 mg/L-250ml-Lithium	Water	250 ml	1000 mg/L
IC500.Li1000.W	Li	1000 mg/L-500ml-Lithium	Water	500 ml	1000 mg/L
IC125.MG1000.W	Mg	1000 mg/L-125ml-Magnesium	Water	125 ml	1000 mg/L
IC250.MG1000.W	Mg	1000 mg/L-250ml-Magnesium	Water	250 ml	1000 mg/L
IC500.MG1000.W	Mg	1000 mg/L-500ml-Magnesium	Water	500 ml	1000 mg/L
IC125.CA1000.W	Ca	1000 mg/L-125ml-Calcium	Water	125 ml	1000 mg/L
IC250.CA1000.W	Ca	1000 mg/L-250ml-Calcium	Water	250 ml	1000 mg/L
IC500.CA1000.W	Ca	1000 mg/L-500ml-Calcium	Water	500 ml	1000 mg/L
IC125.CL1000.W	Cl	1000 mg/L-125ml-Chloride	Water	125 ml	1000 mg/L
IC250.CL1000.W	Cl	1000 mg/L -250ml- Chloride	Water	250 ml	1000 mg/L
IC500.CL1000.W	Cl	1000 mg/L-500ml-Chloride	Water	500 ml	1000 mg/L

Ion Chromatography Single Ion Standards

Ion Chromatography Single Ion Standards

IC125.I1000.W	I	Ion Chromatography 1000mg/L Iodide Standard in Water	Water	125 ml	1000 mg/l
IC500.I1000.W	I	Ion Chromatography 1000mg/L Iodide Standard in Water	Water	500 ml	1000 mg/l
IC125.PO4P1000.W	PO4-P	Ion Chromatography 1000mg/L Phosphate as Phosphorus Standard in Water	Water	125 ml	1000 mg/L
IC250.PO4P1000.W	PO4-P	Ion Chromatography 1000mg/L Phosphate as Phosphorus Standard in Water	Water	250 ml	1000 mg/L
IC500.PO4P1000.W	PO4-P	Ion Chromatography 1000mg/L Phosphate as Phosphorus Standard in Water	Water	500 ml	1000 mg/L
IC125.NH4N1000.W	NH4-N	Ion Chromatography 1000mg/L Ammonium as N Standard in Water	Water	125 ml	1000 mg/L
IC250.NH4N1000.W	NH4-N	Ion Chromatography 1000mg/L Ammonium as N Standard in Water	Water	250 ml	1000 mg/L
IC500.NH4N1000.W	NH4-N	Ion Chromatography 1000mg/L Ammonium as N Standard in Water	Water	500 ml	1000 mg/L
IC125.NO3N1000.W	NO3-N	Ion Chromatography 1000mg/L Nitrate as N Standard in Water	Water	125 ml	1000 mg/L
IC250.NO3N1000.W	NO3-N	Ion Chromatography 1000mg/L Nitrate as N Standard in Water	Water	250 ml	1000 mg/L
IC500.NO3N1000.W	NO3-N	Ion Chromatography 1000mg/L Nitrate as N Standard in Water	Water	500 ml	1000 mg/L
IC500.NO2N1000.W	NO2-N	Ion Chromatography 1000mg/L Nitrite as N Standard in Water	Water	500 ml	1000 mg/L
IC250.NO2N1000.W	NO2-N	Ion Chromatography 1000mg/L Nitrite as N Standard in Water	Water	250 ml	1000 mg/L
IC125.NO2N1000.W	NO2-N	Ion Chromatography 1000mg/L Nitrite as N Standard in Water	Water	125 ml	1000 mg/L
IC125.CN1000.W	CN	Ion Chromatography 1000mg/L Cyanide Standard in Water	Water	125 ml	1000 mg/L
IC250.CN1000.W	CN	Ion Chromatography 1000mg/L Cyanide Standard in Water	Water	250 ml	1000 mg/L
IC500.CN1000.W	CN	Ion Chromatography 1000mg/L Cyanide Standard in Water	Water	500 ml	1000 mg/L

Ion Chromatography Multi-Ion Standards



Ion Chromatography Multi-Ion Standards

Product Code:	IC125.M18002	Matrix:	Water
Chapter:	125ml-Anion Reference Standard Mixture - 18002	Volume:	125 ml
Contents			
F: 1000, Cl: 1000, NO ₂ : 1000, Br: 1000, NO ₃ : 1000, SO ₄ : 1000, PO ₄ : 1000 (mg/l)			

Product Code:	IC125.M18003	Matrix:	Water
Chapter:	125ml-Anion Reference Standard Mixture - 18003	Volume:	125 ml
Contents			
F: 20, Cl: 30, NO ₂ : 100, Br: 100, NO ₃ : 100, SO ₄ : 150, PO ₄ : 150 (mg/l)			

Product Code:	IC125.M18004	Matrix:	Water
Chapter:	125ml-Cation Reference Standard Mixture - 18004	Volume:	125 ml
Contents			
NH ₄ : 100, Ca: 100, Ba: 100, Mn: 100, Sr: 100, Li: 100, Mg: 100, K: 100, Na: 100 (mg/l)			

Product Code:	IC125.M18005	Matrix:	Water
Chapter:	125ml-Cation Reference Standard Mixture 18005	Volume:	125 ml
Contents			
NH ₄ : 400, Ca: 1000, Li: 50, Mg: 200, K: 200, Na: 200 (mg/l)			



The certified values are obtained by gravimetric and volumetric preparations. It is confirmed against certified reference materials traceable to SI of NIST.

Ion Chromatography Multi-Ion Standards

Ion Chromatography Multi-Ion Standards

Product Code:	IC500.M18002	Matrix:	Water
Chapter:	500ml-Anion Reference Standard Mixture - 18002	Volume:	500 ml
Contents			
F: 1000, Cl: 1000, NO ₂ : 1000, Br: 1000, NO ₃ : 1000, SO ₄ : 1000, PO ₄ : 1000 (mg/l)			

Product Code:	IC500.M18003	Matrix:	Water
Chapter:	500ml-Anion Reference Standard Mixture - 18003	Volume:	500 ml
Contents			
F: 20, Cl: 30, NO ₂ : 100, Br: 100, NO ₃ : 100, SO ₄ : 150, PO ₄ : 150 (mg/l)			

Product Code:	IC500.M18004	Matrix:	Water
Chapter:	500ml-Cation Reference Standard Mixture - 18004	Volume:	500 ml
Contents			
NH ₄ : 100, Ca: 100, Ba: 100, Mn: 100, Sr: 100, Li: 100, Mg: 100, K: 100, Na: 100 (mg/l)			

Product Code:	IC500.M18005	Matrix:	Water
Chapter:	500ml-Cation Reference Standard Mixture 18005	Volume:	500 ml
Contents			
NH ₄ : 400, Ca: 1000, Li: 50, Mg: 200, K: 200, Na: 200 (mg/l)			

Eluent Concentrates

Product Code	Symbol	Description	Matrix	Volume	Value
EC.NA ₂ CO ₃ .0.5M		Volumetric Solution for Chromatography 0.5M Na ₂ CO ₃	Water	500 ml	0.5 M
EC.NA ₂ CO ₃ .0.05M		Volumetric Solution for Chromatography 0.05M Na ₂ CO ₃	Water	500 ml	0.05 M
EC.NAHCO ₃ .0.5M		Volumetric Solution for Chromatography 0.5M NaHCO ₃	Water	500 ml	0.5 M
EC.NAHCO ₃ .0.05M		Volumetric Solution for Chromatography 0.05M NaHCO ₃	Water	500 ml	0.05 M
EC.MSA.0.5M		Volumetric Solution for Chromatography 0.5M MSA	Water	500 ml	0.5 M

ANALYTICAL REAGENTS



Analytical Reagents and Solutions

Analytical Reagents and Solutions

pH Standards

Product Code	Symbol	Description	Matrix	Volume	Value
PH500.2.W		pH 2 Reference Standard Solution	Water	500 ml	pH 2
PH1000.2.W		pH 2 Reference Standard Solution	Water	1000 ml	pH 2
PH500.3.W		pH 3 Reference Standard Solution	Water	500 ml	pH 3
PH1000.3.W		pH 3 Reference Standard Solution	Water	1000 ml	pH 3
PH500.4.W		pH 4 Reference Standard Solution	Water	500 ml	pH 4
PH1000.4.W		pH 4 Reference Standard Solution	Water	1000 ml	pH 4
PH500.5.W		pH 5 Reference Standard Solution	Water	500 ml	pH 5
PH1000.5.W		pH 5 Reference Standard Solution	Water	1000 ml	pH 5
PH500.6.W		pH 6 Reference Standard Solution	Water	500 ml	pH 6
PH1000.6.W		pH 6 Reference Standard Solution	Water	1000 ml	pH 6
PH500.7.W		pH 7 Reference Standard Solution	Water	500 ml	pH 7
PH1000.7.W		pH 7 Reference Standard Solution	Water	1000 ml	pH 7
PH500.8.W		pH 8 Reference Standard Solution	Water	500 ml	pH 8
PH1000.8.W		pH 8 Reference Standard Solution	Water	1000 ml	pH 8
PH500.9.W		pH 9 Reference Standard Solution	Water	500 ml	pH 9
PH1000.9.W		pH 9 Reference Standard Solution	Water	1000 ml	pH 9
PH500.10.W		pH 10 Reference Standard Solution	Water	500 ml	pH 10
PH1000.10.W		pH 10 Reference Standard Solution	Water	1000 ml	pH 10
PH500.11.W		pH 11 Reference Standard Solution	Water	500 ml	pH 11
PH1000.11.W		pH 11 Reference Standard Solution	Water	1000 ml	pH 11
PH500.12.W		pH 12 Reference Standard Solution	Water	500 ml	pH 12
PH1000.12.W		pH 12 Reference Standard Solution	Water	1000 ml	pH 12



The reference standard solution is manufactured using a high-purity starting material. The matrix of standard solution is 18 MOhm deionized water and chromatography grade acid solution

Analytical Reagents and Solutions

Analytical Reagents and Solutions

Conductivity Standards

Product Code	Symbol	Description	Matrix	Volume	Value
CY500.500.W		500 uS/cm Conductivity Reference Standard Solution	Water	500 ml	500 uS/cm
CY1000.500.W		500 uS/cm Conductivity Reference Standard Solution	Water	1000 ml	500 uS/cm
CY500.1000.W		1000 uS/cm Conductivity Reference Standard Solution	Water	500 ml	1000 uS/cm
CY1000.1000.W		1000 uS/cm Conductivity Reference Standard Solution	Water	1000 ml	1000 uS/cm
CY500.147.W		147 uS/cm Conductivity Reference Standard Solution	Water	500 ml	147 uS/cm
CY1000.147.W		147 uS/cm Conductivity Reference Standard Solution	Water	1000 ml	147 uS/cm
CY500.1413.W		1413 uS/cm Conductivity Reference Standard Solution	Water	500 ml	1413 uS/cm
CY1000.1413.W		1413 uS/cm Conductivity Reference Standard Solution	Water	1000 ml	1413 uS/cm
CY500.2000.W		2000 uS/cm Conductivity Reference Standard Solution	Water	500 ml	2000 uS/cm
CY1000.2000.W		2000 uS/cm Conductivity Reference Standard Solution	Water	1000 ml	2000 uS/cm
CY500.1500.W		1500 uS/cm Conductivity Reference Standard Solution	Water	500 ml	1500 uS/cm
CY1000.1500.W		1500 uS/cm Conductivity Reference Standard Solution	Water	1000 ml	1500 uS/cm
CY500.3000.W		3000 uS/cm Conductivity Reference Standard Solution	Water	500 ml	3000 uS/cm
CY1000.3000.W		3000 uS/cm Conductivity Reference Standard Solution	Water	1000 ml	3000 uS/cm

Turbidity Standards

Product Code	Symbol	Description	Matrix	Volume	Value
TY500.10.W		10 NTU Turbidity Reference Standard Solution	Water	500 ml	10 NTU
TY100.10.W		10 NTU Turbidity Reference Standard Solution	Water	100 ml	10 NTU
TY500.20.W		20 NTU Turbidity Reference Standard Solution	Water	500 ml	20 NTU
TY100.20.W		20 NTU Turbidity Reference Standard Solution	Water	100 ml	20 NTU
TY500.50.W		50 NTU Turbidity Reference Standard Solution	Water	500 ml	50 NTU
TY100.50.W		50 NTU Turbidity Reference Standard Solution	Water	100 ml	50 NTU
TY100.100.W		100 NTU Turbidity Reference Standard Solution	Water	100 ml	100 NTU
TY500.100.W		100 NTU Turbidity Reference Standard Solution	Water	500 ml	100 NTU
TY500.200.W		200 NTU Turbidity Reference Standard Solution	Water	500 ml	200 NTU
TY100.200.W		200 NTU Turbidity Reference Standard Solution	Water	100 ml	200 NTU
TY500.500.W		500 NTU Turbidity Reference Standard Solution	Water	500 ml	500 NTU
TY100.500.W		500 NTU Turbidity Reference Standard Solution	Water	100 ml	500 NTU
TY500.1000.W		1000 NTU Turbidity Reference Standard Solution	Water	500 ml	1000 NTU
TY100.1000.W		1000 NTU Turbidity Reference Standard Solution	Water	100 ml	1000 NTU
TY500.2000.W		2000 NTU Turbidity Reference Standard Solution	Water	500 ml	2000 NTU
TY100.2000.W		2000 NTU Turbidity Reference Standard Solution	Water	100 ml	2000 NTU

Analytical Reagents and Solutions

Analytical Reagents and Solutions

TOC Standards

Product Code	Symbol	Description	Matrix	Volume	Value
TOC125.10.W		10 mg/l - 125 ml - TOC Standardas potassium hydrogen phthalate	Water	125 ml	10 mg/l
TOC500.10.W		10 mg/l - 500 ml - TOC Standardas potassium hydrogen phthalate	Water	500 ml	10 mg/l
TOC125.100.W		100 mg/l - 125 ml - TOC Standardas potassium hydrogen phthalate	Water	125 ml	100 mg/l
TOC500.100.W		100 mg/l - 500 ml - TOC Standardas potassium hydrogen phthalate	Water	500 ml	100 mg/l
TOC125.200.W		200 mg/l - 125 ml - TOC Standardas potassium hydrogen phthalate	Water	125 ml	200 mg/l
TOC500.200.W		200 mg/l - 500 ml - TOC Standardas potassium hydrogen phthalate	Water	500 ml	200 mg/l
TOC500.500.W		500 mg/l - 500 ml - TOC Standardas potassium hydrogen phthalate	Water	500 ml	500 mg/l
TOC125.500.W		500 mg/l - 125 ml - TOC Standardas potassium hydrogen phthalate	Water	125 ml	500 mg/l
TOC500.1000.W		1000 mg/l - 500 ml - TOC Standardas potassium hydrogen phthalate	Water	500 ml	1000 mg/l
TOC125.1000.W		1000 mg/l - 125 ml - TOC Standardas potassium hydrogen phthalate	Water	125 ml	1000 mg/l

General Reference Standards

General Reference Standards

Miscellaneous Standards

Product Code	Symbol	Description	Matrix	Volume	Value
Setsol.1000		1000 mg/l - Total Settleable Solids Reference Standard	Water	1000 ml	1000 mg/l
TSS.1000		1000 mg/l - Total Suspended Solids Reference Standard.	Water	1000 ml	1000 mg/l
TSM.1000		1000 mg/l - Total Solid Matters Reference Standard	Water	1000 ml	1000 mg/l
SAL.10000		10.000 mg/l - Salinity Standard	Water	500 ml	10.000 mg/l
COD125.1000.W		1000 mg/l - 125 ml - COD Reference Standard (From KHP)	Water	125 ml	1000 mg/l
COD125.500.W		500 mg/l 125 ml COD Reference Standard (From KHP)	Water	125 ml	500 mg/l
COD125.2000.W		2000 mg/l - 125 ml - COD Reference Standard (From KHP)	Water	125 ml	2000 mg/l
COD30x4.500.W		500 mg/l - 30 ml x 4 - COD Reference Standard (From KHP)	Water	30ml x 4	500 mg/l
PL100.1000.W		Phenol 1000 mg/l Reference Standard Solution in Water	Water	100 ml	1000 mg/l
PL500.1000.W		Phenol 1000 mg/l Reference Standard Solution in Water	Water	500 ml	1000 mg/l
MBASA100.1000.W		MBAS-Alkyl 1000 mg/l Reference Standard Solution in Water	Water	100 ml	1000 mg/l
MBAS100.1000.W		MBAS-Alkylbenzene 1000 mg/l Reference Standard Solution in Water	Water	100 ml	1000 mg/l
MBAS30x4.1000.W		MBAS-Alkylbenzene 1000 mg/l Reference Standard Solution in Water	Water	30ml x 4	1000 mg/l
MBASA30x4.1000.W		MBAS-Alkyl 1000 mg/l Reference Standard Solution in Water	Water	30ml x 4	1000 mg/l
TDS.1000		1000 mg/l - Total Dissolved Solids Reference Standard	Water	1000 ml	1000 mg/l
FCL100.1000.W		1000 mg/l - Total Residual Free Chlorine Reference Standard Solution in Water	Water	100 ml	1000 mg/l
FCL500.1000.W		1000 mg/l - Total Residual Free Chlorine Reference Standard Solution in Water	Water	500 ml	1000 mg/l
OG25.200.ACE		200 mg/l - 25 ml - Oil&Grease Reference Standard Solution in Acetone (Certified value is obtained after diluton to 1 liter water)	Acetone	25 ml	200 mg/l
OG100.50.WtrACE		50 mg/l - 100 ml - Oil&Grease Reference Standard Solution in Water (Ready to Use)	Water	100 ml	50 mg/l
TCN100.1000.W		1000 mg/l - 100 ml -Total Cyanide Reference Standard Solution in Water	Water	100 ml	1000 mg/l
TCN500.1000.W		1000 mg/l - 500 ml -Total Cyanide Reference Standard Solution in Water	Water	500 ml	1000 mg/l
TN100.1000.W		1000 mg/l - 100 ml - Total Nitrogen Reference Standard Solution in Water	Water	100 ml	1000 mg/l
TN500.1000.W		1000 mg/l - 500 ml - Total Nitrogen Reference Standard Solution in Water	Water	500 ml	1000 mg/l
TKN100.1000.W		1000 mg/l - 100 ml - Total Kjeldahl Nitrogen Reference Standard Solution in Water	Water	100 ml	1000 mg/l
TKN500.1000.W		1000 mg/l - 500 ml - Total Kjeldahl Nitrogen Reference Standard Solution in Water	Water	500 ml	1000 mg/l
TP100.1000.W		1000 mg/l - 100 ml - Total Phosphorus Reference Standard Solution in Water	Water	100 ml	1000 mg/l
TP500.1000.W		1000 mg/l - 500 ml - Total Phosphorus Reference Standard Solution in Water	Water	500 ml	1000 mg/l
TPL100.1000.W		1000 mg/l - 100 ml - Total Phenolics Reference Standard Solution in Water	Water	100 ml	1000 mg/l
TPL500.1000.W		1000 mg/l - 500 ml - Total Phenolics Reference Standard Solution in Water	Water	500 ml	1000 mg/l
TCR100.1000.W		1000 mg/l - 100 ml - Total Cyanide Reference Standard Solution in Water	Water	100 ml	1000 mg/l

General Reference Standards

General Reference Standards

TCR500.1000.W		1000 mg/l - 500 ml - Total Cyanide Reference Standard Solution in Water	Water	500 ml	1000 mg/l
AK100.1000.W		1000 mg/l - 100 ml - Alkalinity Reference Standard Solution (As CaCO ₃) in Water	Water	100 ml	1000 mg/l
AK500.1000.W		1000 mg/l - 500 ml - Alkalinity Reference Standard Solution (As CaCO ₃) in Water	Water	500 ml	1000 mg/l
AK500.1000.W		1000 mg/l - 30 ml x 4 - Calcium Hardness Reference Standard Solution in Water	Water	30ml x 4	1000 mg/l
HS30x4.1000.1CA		1000 mg/l - 100 ml - Calcium Hardness Reference Standard Solution in Water	Water	100 ml	1000 mg/l



All raw materials used for manufacture of reference materials are selected at high purity level. Following completion of manufacture, bottles are made ready for shipment by being placed in aluminum locked packages to minimize environmental effects and for extra protection.

STANDARDS FOR LEGAL REGULATIONS



Calibration Standards for Legal Regulations

Calibration Standards for Legal Regulations

Product Code:	ICP250.SET1	Matrix:	Water
Chapter:	IC Anion Calibration Reference Set - 1	Volume:	250 ml
Contents			
Water analysis reference standard set for 7 point calibration curve. Calibration set contains 250 ml x 7 bottles.			
1-) 250 ml / F - 0,25 mg/l; PO4 - 0,3 mg/L; NO3 - 5 mg/L; Cl - 10 mg/L; SO4 - 10 mg/L			
2-) 250 ml / F - 0,5 mg/l; PO4 - 0,5 mg/L; NO3 - 10 mg/L; Cl - 20 mg/L; SO4 - 20 mg/L			
3-) 250 ml / F - 0,75 mg/l; PO4 - 0,75 mg/L; NO3 - 20 mg/L; Cl - 30 mg/L; SO4 - 30 mg/L			
4-) 250 ml / F - 1 mg/l; PO4 - 1 mg/L; NO3 - 25 mg/L; Cl - 40 mg/L; SO4 - 40 mg/L			
5-) 250 ml / F - 1,25 mg/l; PO4 - 1,5 mg/L; NO3 - 35 mg/L; Cl - 50 mg/L; SO4 - 50 mg/L			
6-) 250 ml / F - 1,50 mg/l; PO4 - 2,0 mg/L; NO3 - 40 mg/L; Cl - 60 mg/L; SO4 - 60 mg/L			
7-) 250 ml / F - 2,0 mg/l; PO4 - 2,5 mg/L; NO3 - 45 mg/L; Cl - 70 mg/L; SO4 - 70 mg/L			

Product Code:	ICP250.SET2	Matrix:	Water
Chapter:	IC Cation Calibration Reference Set - 2	Volume:	250 ml
Contents			
Water analysis reference standard set for 6 point calibration curve. Calibration set contains 250 ml x 6 bottles.			
1-) 250 ml / K - 1 mg/L; Ca - 5 mg/L; Mg - 5 mg/L; Na - 5 mg/L			
2-) 250 ml / K - 3 mg/L; Ca - 10 mg/L; Mg - 10 mg/L; Na - 10 mg/L			
3-) 250 ml / K - 5 mg/L; Ca - 20 mg/L; Mg - 15 mg/L; Na - 20 mg/L			
4-) 250 ml / K - 10 mg/L; Ca - 30 mg/L; Mg - 20 mg/L; Na - 30 mg/L			
5-) 250 ml / K - 15 mg/L; Ca - 40 mg/L; Mg - 30 mg/L; Na - 40 mg/L			
6-) 250 ml / K - 20 mg/L; Ca - 45 mg/L; Mg - 35 mg/L; Na - 45 mg/L			

CROSS REFERENCE STANDARDS



Cross Reference Standards

Cross Reference Standards

Product Code:	ICP125.M28000	Matrix:	10% HNO ₃
Chapter:	Equivalent to Agilent	Volume:	125 ml
Contents			
Equivalent Code: 5188-6525 mg/l Bi: 10, Ge: 10, In: 10, Li-6: 10, Lu: 10, Rh: 10, Sc: 10, Tb:10			


Product Code:	ICP125.M28015	Matrix:	2% HNO ₃
Chapter:	Equivalent to Agilent	Volume:	125 ml
Contents			
Equivalent Code: 5185-5959 ug/l Ce: 1, Co: 1, Li: 1, Mg: 1, Tl: 1, Y: 1			

Product Code:	IC125.M18000	Matrix:	Water
Chapter:	Equivalent to Thermo / Dionex	Volume:	125 ml
Contents			
Ion Chromatography 7 anion standard mg/l F: 20, Cl: 100, NO ₂ : 100, Br: 100, NO ₃ : 100, SO ₄ : 100, PO ₄ : 200			

Product Code:	IC125.M18001	Matrix:	Water
Chapter:	Equivalent to Thermo / Dionex	Volume:	125 ml
Contents			
Ion Chromatography 6 cation Standard mg/l Li: 50, Na: 200, NH ₄ : 250, K: 500, Ca: 500, Mg: 250			

CERTIFICATION PROCESS

- Purity of starting material shown in the certificates.
- Preparation and certification informations are enclosed.
- Traceability information is stated.
- Results of stability and homogeneity studies are noted.
- Expiration of the product is shown.
- Barcode of solution is identified in the certificate.



Address: Mahmutiye M. Bahadır C. No:27 / 34 Mirsin - Turkey
Phone: +90850 303 40 76 / Fax: +90850 203 51 53
E-Mail: global@labsert.com / Web: www.labsertchemical.com

CERTIFICATE OF ANALYSIS

1. COMPANY INFORMATIONS

Labsert is an ISO 9001 registered company which manufactures reference materials for international markets. All tests and analysis are carried out by ISO/IEC 17025 "General Requirements for the Competence of Testing and Calibration Laboratories" accredited laboratory as a Labsert's sub-contracted business partner.

2. DESCRIPTION OF CERTIFIED REFERENCE MATERIAL (CRM)

Product Name	ICP-MS - Mn (Manganese) Reference Standard Solution		
Catalog Code	ICP-Mn1000_2NA		
Lot Number	ICP-Mn1000- xxx		
Concentration	1000 mg/l		
Starting Material	NH ₄ H ₂ PO ₄		
Starting Material Purity* (%)	99,90		
Matrix	2% HNO ₃		
Density (g/cm ³ at 25 °C)	0,9998		

3. CERTIFIED VALUE / UNCERTAINTY

CAS	Component	Symbol	Certified Value	Uncertainty ±	Calculation
7439-96-5	Manganese	Mn	998 mg/l	2,8 mg/l	Mean of the Results

4. PREPARATION AND CERTIFICATION INFORMATIONS

4.1. The certified value is obtained by gravimetric and volumetric preparations. It is confirmed against certified reference materials traceable to SI of NIST (National Institute of Standards and Technology), using one or two instrumental methods by ISO 17025 accredited laboratory which is sub-contracted business partner of Labsert. Certified value where one method is used, is the mean of individual results. Certified value where two methods are used, is the weighted mean of two results. The reported expanded uncertainty of measurement is calculated for a 95% confidence interval using a coverage factor of $k = 2$.

4.2. The reference standard solution is manufactured using a high-purity starting material. The matrix of standard solution is 18 MOhm deionized water and chromatography grade acid solution (if necessary).

5. TRACEABILITY STATEMENT

- Property of the result of a measurement or the value of a standard whereby it can be related to stated references, usually national or international standards, through an unbroken chain of comparisons all having stated uncertainties (ISO VIM, 2nd ed., 1993, definition 6.10)
- The measurement results are traceable to SI. All analytical balances are calibrated yearly by ISO 17025 accredited metrology service and are verified monthly by an in-house method using NIST traceable analytical weights. Analytical balances are also weekly checked by using Class A laboratory glassware.
- Automatic equipments are calibrated yearly verified by ISO 17025 accredited metrology service and are verified weekly by an in-house method.
- The results from temperature measurement are traceable to SI. Thermometers are calibrated yearly by ISO 17025 accredited metrology service.

6. INTENDED USE

This reference standard solution is for laboratory use only and it is intended for the calibration of analytical instruments such as Ion Chromatography, HPLC, UV/VIS, MS, ICP. This reference standard solution is also intended for validation of analytical methods, detection limit and linearity studies, preparation of "working reference samples". This statement is not intended to restrict the use for other purposes.

* These values are not certified

7. INSTRUCTIONS FOR THE CORRECT USE OF THIS REFERENCE MATERIAL

7.1. Stability, Storage and Handling

This certified reference material is used with a guaranteed stability until $\pm 0.5\%$ of the certified concentration within its shelf life. For the guaranteed stability, keep it in its original packaging, tightly capped when it is not in use and store at $20 \pm 4^\circ\text{C}$. Do not pipette from the container. Do not return portions removed from pipetting to container.

8. HAZARDOUS INFORMATION

Please refer to Safety Data Sheet (SDS) for hazardous information regarding this certified reference material.

9. HOMOGENEITY

This solution was mixed according to an in-house procedure and is guaranteed to be homogeneous. To achieve the sufficient homogeneity, please mix the sample slowly by inversion while it is capped.

10. EXPIRATION INFORMATIONS

The expiration date is guaranteed to be valid when this certified reference material is kept in original unopened bottle and aluminized bag according to the instructions in article 7.1. Recommended period of use should not exceed 6 months from the date of opening.

Certification Date: xxx
Retest Date:
Shipment Date (MM.YY): 04.2018
Expiration Date: 04.2020

11. NAMES OF CERTIFYING OFFICERS

Prepared By:

Controlled By:



TRADING PROGRESS

Quotation and Ordering

Regarding international order progress, you get our quotation first to see the total cost (Including shipment cost, packaging cost etc.) before payment on our web site. To get our quotation:

1. Select the products you are interested in and add them into your shopping cart. On the Check-Out page, please select Bank Transfer option and confirm your order to get our quotation by e-mail which includes all related costs and "payment page" link. (During this progress, we do not request payment from you. This process is just for getting our quotation)
2. When you click "payment page" link in the e-mail, our web page regarding international payments is opened to let you choose the best method between credit card payment and bank transfer payment.
3. Then you can select the best payment method for you on the page. You send the total cost according to our quotation. Once your payment is received by us, purchasing progress begins. We inform you about manufacturing and shipment progress.

You can also send us your request by e-mail

For international orders: global@labsert.com / For domestic orders: bilgi@labsert.com

Conditions of Sales and Shipping

Prices are in US dollars or Turkish Lira and are subject to change without notice. Orders may be placed via the telephone, fax, mail, or web site. All orders must clearly indicate the customer telephone number, fax number, complete billing address, complete shipping address, catalog numbers, product names, quantity ordered and price. Customers are able to pay via credit card or bank transfer on our web site. Payment must be done to start the purchasing process.

For international orders, products are sent FOB shipping point. DHL, TNT UPS, etc., shipping charges are prepaid and added to your invoice. Handling fee may be applied to each order. International Customers select bank transfer for the payment option. Then we send the exact quotation with the all related costs and customers pay according to our quotation via credit card or bank transfer before the shipment.

For domestic orders, products are sent by local cargo companies or your contracted cargo company as cash on delivery.

For more information, please check the link below.

<https://www.labsertchemical.com/Ordering-Delivery-Information>



*International Manufacturer of Reference
Standards and Analytical Reagents*

 /Labsertchemical

 /Labsert

 /Labsert

 /labsert

Labsert Laboratory Services and Technologies

Address: Mahmudiye M. Matiat Is Mrk. No:5 / 34 Akdeniz-Mersin-Turkey
Phone: +908503034076
Fax: +908502035153
Web: www.labsertchemical.com
E-Mail: global@labsert.com