



Questron Technologies Corp.

Innovative Sample Preparation Solutions for Analytical Laboratories

Automated Hot Block Digestion

Vulcan Automated Digestion and Work-up Station

- As many as 84 samples can be automatically digested and worked up
- A true *click-and-walk-away* system – frees up valuable chemists' time to perform other tasks
- Total PC control helps storage of digestion and sample rack preparation data
- PC-controlled preparation of sample racks
- 6 dedicated lines for reagents
- Complete and thorough mixing of samples by stream of air bubbles
- Leveling up after digestions
- Built-in fumehood
- Also available: Exhaust fume scrubbers, exhaust blowers and reagent storage cabinets – all-plastic construction and compact in size to fit in any laboratory



Microwave Digestion Systems

Versatility, Safety and Productivity in Closed Vessel Digestions

- Configurations to suit the needs and budgets of all laboratories
- Multilayered safety features employing our unique OPGuard™ and Q-Sens™ technologies
- Corrosion resistant PTFE coated reinforced cavity
- Spring-loaded door latch providing effective safety shield
- Door safety interlocks
- Vessel sets available for high pressure and high throughput

QLAB Pro

- Ideal for digestion, extraction, evaporation and synthesis
- Built-in temperature monitoring of digestion vessels via IR sensing
- User friendly control panel, operated through full-colour bright LCD touchscreen
- Preloaded with EPA and other standard methods

- Choice of sensor vessels for in-situ temperature and pressure measurement
- Can be used as a stand-alone system or with PC
- On-board storage of 10 time-to-temperature and 5 time-to-power multi-step digestion recipes
- Real-time graphical display for digestion parameters

QWAVE 1000

- Four power level settings
- Sample digestion based on time-to-power control
- Preferred choice for high-throughput laboratories digesting known samples of all sizes
- Practical and cost-effective alternative for conventional digestion techniques

Microwave Digestion Vessels

Ultra High Pressure (UHP), Very High Pressure (VHP) and Low Volume High Pressure (LVHP)

Vessel	Vessels/Carousel	Total Volume/Vessel	Max. Control Press.	Max. Control Temp.
UHP ¹	10	100ml	1500psi	230°C
VHP ²	10	110ml	625psi	230°C
LVHP	30	50ml	300psi	200°C

1. Data represent both standard digestion vessels and sensor vessels.

2. Data represent both standard digestion vessels and sensor vessels.

Microwave Ashing System

QAsh 1800

- Faster, cleaner and reproducible computer-controlled ashing
- Cuts down on ashing times for most samples by 90% or more
- Display and storage of sample ashing data in tabular and graphical format
- Weighing balance integration capability




Hot Block Digestion

QBlock Digestion System – with Quad-channel Heating Controller



QBlock Commander

- Controls up to four QBlock Digesters or third-party hot blocks independently
- Operates from outside the fumehood in acid free environment
- Stores up to 10 multistep digestion programs
- PC hook-up capability for digestion data transfer



QBlock Digesters

- Capable of reaching up to 230°C for higher temperature digestions
- Teflon®-coated graphite block and exterior
- Models available to accommodate 15, 50 or 100ml vials
- Polypropylene, Teflon® and glass vials available



Available Vials

	Polypropylene	Glass	Teflon®
50ml	Yes	Yes	Yes
100ml	Yes	No	Yes
Temp.	120°C	230°C	230°C
Usage	Single	Multiple	Multiple
Cap	Yes	No	No

Your local distributor:

Qestron Technologies Corp.

6725 Millcreek Drive, Unit 7, Mississauga, Ontario L5N 5V3, Canada
Tel: 905-363-1223 Fax: 905-363-1227 E-Mail: info@qtechcorp.com Website: www.qtechcorp.com